## User's Instructions

| Catalogue number | Size                | Scan to documents                                                                   |
|------------------|---------------------|-------------------------------------------------------------------------------------|
| 65 8121          | Extra Small (XS/6)  |  |
| 65 8122          | Small (S/7)         |                                                                                     |
| 65 8123          | Medium (M/8)        |                                                                                     |
| 65 8124          | Large (L/9)         |                                                                                     |
| 65 8125          | Extra Large (XL/10) |                                                                                     |

|                                                                                   |                                                                |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------|
|  | EN ISO 21420:2020 = Prior to use, please read the instructions |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------|

|                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br>2797 | The product complies with and is certified in accordance with the requirements of EU Regulation 2016/425. This product is EU conformity to type based on internal production control plus supervised product checks at random intervals (Module C2) by BSI Group, The Netherlands B.V. (Notified Body No: 2797), Say Building - John M. Keynesplein 9, Unit 4.2, 1066 EP Amsterdam - The Netherlands |
|--------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br>0086 | The product complies with and is certified in accordance with the requirements of Regulation 2016/425 on personal protective equipment, as amended to apply in Great Britain. This product is UKCA conformity to type based on internal production control plus supervised product checks at random intervals (Module C2) by BSI Assurance UK Limited (Approved Body No: 0086) Kitemark Court - Davy Avenue - Knowhill - Milton Keynes - MK5 8PP - United Kingdom |
|---------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

| Determination of material resistance to permeation by chemicals (EN ISO 374-1:2016 + A1:2018) and degradation (EN ISO 374-4:2019)                                       |                          |            |                                              |                            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|------------|----------------------------------------------|----------------------------|
|   | Code letter and chemical | CAS number | Permeation performance and breakthrough time | Degradation (DR) Mean in % |
|                                                                                                                                                                         | K: Sodium hydroxide 40%  | 1310-73-2  | Level 6, > 480 min.                          | -10.8                      |
|                                                                                                                                                                         | P: Hydrogen peroxide 30% | 7722-84-1  | Level 6, > 480 min.                          | -5.9                       |
|                                                                                                                                                                         | T: Formaldehyde 37%      | 50-00-0    | Level 6, > 480 min.                          | -37.7                      |

| Marking of gloves protecting against viruses, bacteria and fungi, as defined in EN ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms – Part 5: Terminology and performance requirements for micro-organisms risks). Protection against viruses: tested according to ISO 16604:2004 (Procedure B) Phi-X, Herpes, HIV. |                   |                               |                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|-------------------------------|------------------|
|                                                                                                                                                                                | Performance level | Acceptable quality level unit | Inspection level |
|                                                                                                                                                                                                                                                                                                                                                      | Level 3           | < 0.65                        | G1               |
|                                                                                                                                                                                                                                                                                                                                                      | Level 2           | < 1.5                         | G1               |
|                                                                                                                                                                                                                                                                                                                                                      | Level 1           | < 4.0                         | S4               |

(en) Storage recommendations: Store in a cool and dry place. Storage temperature 5°C (41°F) to 40°C (104°F). Gloves should be protected from exposure to direct sunlight, intense artificial light, X-ray machines and other sources of ozone. Recommendation for packaging whilst being transported: Gloves should be transported in their original packaging to maintain the protective properties and usability characteristics of the gloves. Disposal: Dispose of in accordance with local authority regulations. This is an MDR Class 1 examination glove according to Regulation (EU) 2017/745, a PPE Category III protective glove according to Regulation (EU) 2016/425 and Regulation 2016/425 on Personal Protective Equipment, as amended to apply in Great Britain. Gloves provide protection against dangerous chemicals and micro-organisms risks. Caution: These gloves should not be used for mechanical and/or thermal protection. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only [except in cases where the glove is equal to or over 400 mm - where the cuff is tested also] and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. DONNING: 1. Clean and prepare the work surface before opening the packaging. 2. Remove first glove from its packaging and check for any possible mechanical damage like holes or tears before use. 3. Touching the inner cuff area, pick up the glove and don pulling it all the way up to the wrist or forearm. Avoid any contact with the outer area of the glove. 4. Ensure that the glove fits properly. 5. Repeat steps 1-4 for second glove. DOFFING: 1. Hold the outer edge of the gloves near the cuff/ wrist or forearm area. 2. Carefully pull the glove away from the hand, turning it inside out. Hold it in the opposite hand. 3. Slip two fingers into the cuff of the remaining gloved hand without touching the outer area of the glove. 4. Peel off the glove, inside out, until it is completely removed ensuring that the glove removed encapsulates the first glove in its palm. 5. Discard both gloves appropriately. Declaration of conformity accessible from: <http://shieldscientific.com/resource-center/glove-declaration-of-conformity>.

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**SHIELD Scientific B.V.** Costerweg 1-B, 6702 AA Wageningen, The Netherlands

| EU                                                                                                                                                                                                              | UK                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>NOTIFIED BODY NO: 2797</b><br><b>BSI Group The Netherlands B.V.</b><br>Say Building, John M. Keynesplein 9, Unit 4.2,<br>1066 EP Amsterdam, The Netherlands                                                  | <b>APPROVED BODY NO: 0086</b><br><b>BSI Assurance UK Ltd</b><br>Kitemark Court, Davy Avenue, Knowlhill, Milton<br>Keynes, MK 5 8PP, United Kingdom                                                  |
| <b>MODULE B EU TYPE EXAMINATION UNDER-<br/>TAKEN BY NOTIFIED BODY NO: 2797</b><br><b>BSI Group The Netherlands B.V.</b><br>Say Building, John M. Keynesplein 9, Unit 4.2,<br>1066 EP Amsterdam, The Netherlands | <b>MODULE B UKCA TYPE EXAMINATION<br/>UNDERTAKEN BY APPROVED BODY NO: 0086</b><br><b>BSI Assurance UK Ltd</b><br>Kitemark Court, Davy Avenue, Knowlhill, Milton<br>Keynes, MK 5 8PP, United Kingdom |

**“duoSHIELD™, for hand that make a difference”**

[www.shieldscientific.com](http://www.shieldscientific.com)