

HighPrep™ DTR

Dye Terminator Removal - Bead Based Sequencing Clean Up System

Hassle free replacement of competitor's product for just a fraction of the cost.

Description

HighPrep™ DTR is a paramagnetic bead-based system, designed to remove unincorporated terminators from Sanger sequencing reaction. The HighPrep™ DTR system consists of a selective binding of DNA to the HighPrep™ DTR particles, followed with washing off nucleotides, primers and non-targeted amplicons, and finally elution of pure DNA to be directly used in downstream application. HighPrep™ DTR is designed for both manual and fully automated purification of sequencing products.

Overall Benefits

- Cost Effective. Save up to 30% to similar competitor.
- Get equal or better read lengths than competitor. See Figure 1.
- Long read lengths: Phred20 scores over 700
- Get pass rates of 85% or higher
- No centrifugation or vacuum steps
- Compatible with BigDye® v1.0, 1.1, 2.0, 3.0, 3.1 and other commonly used dyes
- Reduce BigDye® dilution in 1/4, 1/8, 1/16 and 1/32. Reduce sample cost.
- Flexible - can be adapted to most standard liquid handling robot
- Can be used in 96 or 384 well format



Comparative Analysis of HighPrep™ DTR vs CleanSEQ®

Plasmid DNA*	Signal Strength			Reaction Template	
	HighPrep DTR	CleanSEQ	CRL	QV	Ratio
	2155	2096	636	56	0.9375
			635	55	0.125
PCR: 1.8KB PCR fragments*	Signal Strength			Reaction Template	
	HighPrep DTR	CleanSEQ	CRL	QV	Ratio
	1044	872	844	50	1.9375
			776	50	1

Figure 1.
*All samples (plasmid and PCR products) are a median of 16 samples. 5 µl reactions.
NB: No dye blob or ski stop effect were observed.

Ordering Information

Cat No.	Product	Preps
220021	HighPrep DTR (5 mL)	500 ¹ 1,000 ²
220022	HighPrep DTR (50 mL)	5,000 ¹ 10,000 ²
220030	HighPrep DTR (500 mL)	50,000 ¹ 100,000 ²

¹Based on typical 10 µl reaction volume in a 96 well format.

²Based on typical 5 µl reaction volume in a 384 well format.



Biozym Scientific GmbH – Steinbrinksweg 27 - D-31840 Hessisch Oldendorf
Tel.: +49 (0)5152-9020 - support@biozym.com
Biozym Biotech Trading GmbH – Wehlist. 27b – A-1200 Wien
Tel.: +43 (1)3340156-0 - wien@biozym.com
www.biozym.com



USA/Canada: +1 855 262-4246 / +1 301-302-0144

www.magbiogenomics.com

© 2012 MagBio Genomics, Inc. All rights reserved. HighPrep DTR is a trademark of MagBio Genomics, Inc. BigDye is a registered trademark of Applied Biosystems, Inc. /Life Technologies. CleanSEQ is a registered trademark of Beckman Coulter. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.