

# Biozym Probe qPCR Mix Separate ROX

## Product information

Store all kit components at -20°C with minimal exposure to light.

The kit may be stored at 4°C for short term use (up to 1 month). Multiple freeze/thaw cycles (up to 30) are possible with no loss of activity.

Component	331455X 10 rxn	331455S 100 rxn	331455L 500 rxn	331455XL 2500 rxn
Probe qPCR 2x Mix	1 x 0.2 ml	1 x 1.0 ml	5 x 1.0 ml	25 x 1.0 ml
50 µM ROX Additive	1 x 0.05 ml	1 x 0.2 ml	5 x 0.2 ml	25 x 0.2 ml

For research use only.

## 1. Introduction

Biozym Probe qPCR Mix is designed for quantitative, real-time analysis of DNA and cDNA samples. Used polymerase technology along with progressive buffer chemistry enable fast protocols, enhanced specificity and reaction efficiency.

Biozym Probe Mixes are designed for use with probe technologies including TaqMan®, Scorpions® and molecular beacons. The mix is the ideal choice for multiplex qPCR applications.

The used hot-start technology inhibits DNA polymerase activity at ambient temperature, thus preventing formation of primer-dimers to improve reaction sensitivity and specificity.

The reaction chemistry is applicable to most qPCR instruments. ROX Additive is added separately if the instrument requires such as passive reference dye.

## 2. Notes

- Master Mix ideal for multiplex qPCR applications.
- ROX passive reference dye is included separately (see 5.3).

## 3. Technical support

For technical support please contact [support@biozym.com](mailto:support@biozym.com)

### Biozym Scientific GmbH

Steinbrinksweg 27  
D-31840 Hess. Oldendorf  
Ph: +49 5152 9020  
[www.biozym.com](http://www.biozym.com)

#### 4. Reaction conditions for qPCR

Carefully mix and centrifuge the tubes before opening to ensure homogeneity and improve recovery.

**Table 1. Pipetting instructions**

Component	20 µl reaction	Final concentration	Notes
2x qPCR Probe Mix	10 µl	1x	
Forward primer (10 µM)	0.8 µl	400 nM	See 5.2.
Reverse primer (10 µM)	0.8 µl	400 nM	
Probe (10 µM)	0.4 µl	200 nM	See 5.3.
Template DNA cDNA gDNA	<100 ng <1 µg	Variable	
PCR grade water	Up to 20 µl reaction volume		

**Table 2. Cycling instructions**

Cycles	Temperature	Time	Notes
1	95°C	2 min	Initial denaturation and enzyme activation, 2 min for cDNA, 3 min for gDNA
40	95°C 60 to 65°C	5 seconds 20 to 30 seconds	Denaturation Annealing/Extension, do not exceed 30 seconds or temperatures below 60°C

#### 5. Notes about reaction components and cycling conditions

##### 5.1. Biozym Probe qPCR Mix

The 2x mix contains a hot-start DNA polymerase, buffer and dNTPs.

##### 5.2. Primers

For best efficiency of the reaction the amplicon length should be between 80 and 200 bp. Short amplicons enable fast cycling conditions. Primers should have a calculated melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>). The final primer concentration in the reaction should be 0.4 µM (0.2 – 0.4 µM).

##### 5.3. ROX passive reference dye

###### Instrument compatibility

Different real-time PCR instruments may require different levels of ROX passive reference.

###### Addition of ROX additive

The 50 µM ROX Additive supplied is formulated to be added directly to the 1 ml tube of Biozym master mix supplied. Once the ROX is added, the reagent may be directly used or stored at -20°C. Please follow the below table to add the correct amount of ROX for your instrument. Mix carefully after ROX addition.

**Table 3. ROX for Hi-ROX instruments**

Reagent	Hi-ROX instruments	Final concentration	Reaction concentration
Probe qPCR 2x Mix	1.0 ml (0.2 ml)	2x	1x
50 µM ROX Additive	35.0 µl (7 µl)	1.75 µM	875 nM

**Table 4. ROX for Lo-ROX instruments**

Reagent	Lo-ROX instruments	Final concentration	Reaction concentration
Probe qPCR 2x Mix	1.0 ml (0.2 ml)	2x	1x
50 µM ROX Additive	4.0 µl (0.8 µl)	200 nM	100 nM