


riboPOOLs for FFPE samples

Efficient Ribosomal RNA Depletion For FFPE samples and degraded RNA samples


- ✓ Gapless coverage of rRNA
- ✓ Best solution for low-quality RNA (e.g. FFPE)
- ✓ Enzyme-free, fast & easy workflow
- ✓ Hybridization based rRNA removal

Total RNA extracted from FFPE tissue (Formalin-fixed, paraffin-embedded tissues) are heavily degraded with RIN values below 2. Due to their strong fragmentation ribosomal RNAs are difficult to remove before RNA sequencing (RNA-Seq). With a gapless probe coverage of the entire rRNA sequence, siTOOLs Biotech developed an efficient rRNA depletion tool for FFPE and other strongly degraded RNA samples called **FFPE riboPOOL**.



Representative Bioanalyzer Analysis of RNA with varying RNA quality.

RNA-Seq Workflow:



Application:

FFPE riboPOOLs are suitable for rRNA depletion in RNA samples with varying RNA quality (i.e. FFPE and intact RNA samples).

Currently available:



Human FFPE riboPOOL

ID: 57

coverage of entire human rRNA sequence for even rRNA depletion in FFPE and degraded RNA samples



Mouse/Rat FFPE riboPOOL

ID: 58

coverage of entire mouse/rat rRNA sequence for even rRNA depletion in FFPE and degraded RNA samples



Human/Mouse/Rat FFPE riboPOOL

ID: 56

equimolar mixture of Human FFPE riboPOOL & Mouse/Rat FFPE riboPOOL

rRNA depletion efficiency:

The FFPE riboPOOL have been tested on degraded RNA samples with rRNA depletion efficiencies ranging between 87-98%, while the integrity of the transcriptome was maintained (Spearman correlation ≥ 94 between non-depleted and depleted samples).

Further information:

- Enzyme-free depletion
- rRNA depletion completed in 70 mins
- Suitable for RNA input amounts from 10 ng - 3 μ g
- Compatible with various RNA clean-up methods & RNA-Seq library preparation kits

Available formats:

1. Probes alone with nuclease-free water	2 nmol (20 rx)	5 nmol (50 rx)	10 nmol (100 rx)
2. Complete riboPOOL kits (includes buffers, streptavidin-magnetic beads, reaction tubes and ethanol precipitation reagents)	6 rx Trial (one-time purchase)	12 rx	24 rx
			96 rx

Contact us or our [Distributors](#) for questions and orders.
Quick order now on [webshop!](#)

Biozym
SCIENCE IS OUR BUSINESS

Biozym Scientific GmbH
Tel.: 05152/9020, Fax: 05152/2070
Mail: support@biozym.com



www.biozym.com

Biozym Biotech Trading GmbH
Tel.: 01/33401560, Fax: 01/334015688
Mail: wien@biozym.com

