



**nexttec**<sup>TM</sup>  
DNA isolation systems

## Protocol

**DNA Isolation  
from Blood (200 µl)  
by nexttec<sup>TM</sup> 1-Step**

- nexttec<sup>TM</sup> cleanPlate96 -

|                         |                |
|-------------------------|----------------|
| <b>Cat. No. 50N.901</b> | <b>395091N</b> |
| <b>Cat. No. 50N.902</b> | <b>395092N</b> |
| <b>Cat. No. 50N.904</b> | <b>395094N</b> |
| <b>Cat. No. 50N.924</b> | <b>395096N</b> |

Version 3.0

**For research only**

**Biozym**  
SCIENCE IS OUR BUSINESS

**Principle**

nexttec™ 1-Step is the easiest handling and fastest DNA purification system containing a single buffer system and a one-step DNA purification after lysis.

Proteins, detergents and low molecular weight compounds are retained by the nexttec™ sorbent. DNA passes through the nexttec™ cleanColumn during a short, one-step purification procedure.

The obtained DNA is suitable for all common enzymatic reactions (restriction digests, real-time PCR, PCR, genotyping etc.).

**Kit contents**

The kit contains all necessary reagents for lysis and subsequent DNA purification.

| Component                | Art.No.<br>50N.901 | Art.No.<br>50N.902 | Art.No.<br>50N.904 | Art.No.<br>50N.924 |
|--------------------------|--------------------|--------------------|--------------------|--------------------|
| Solution B               | 2 x 120 ml         | 2 x 250 ml         | 2 x 420 ml         | 10 x 500 ml        |
| Proteinase K             | 2.2 ml             | 6 ml               | 9 ml               | 55 ml              |
| Prep Solution            | 40 ml              | 100 ml             | 150 ml             | 2 x 450 ml         |
| DTT                      | 0.2 ml             | 0.5 ml             | 1.0 ml             | 5 ml               |
| Buffer R                 | 14 ml              | 38                 | 57 ml              | 344 ml             |
| nexttec™ cleanPlates96   | 1                  | 2                  | 4                  | 24                 |
| nexttec™ deep-well plate | 3                  | 6                  | 12                 | 72                 |
| Sealing tapes            | 2                  | 4                  | 8                  | 48                 |
| Alu sealing tapes        | 2                  | 4                  | 8                  | 48                 |

**nexttec™ service**

To extend the application range to samples, which are difficult to lyse by the standard procedure, it is recommended to include optional components in the lysis buffer and to optimize the lysis conditions. Please get in contact with [service@nexttec.biz](mailto:service@nexttec.biz) for detailed information.

**Storage Conditions**

During shipment all kit components are stable at room temperature. After arrival, store **Proteinase K, Solution B, Buffer R** and **Prep Solution** at **+2 °C to +8 °C**. Store **DTT** after first opening at **-18 °C to -25 °C**.

**nexttec™ cleanPlates96** are stored at **room temperature (+20 °C to +25 °C)**. If properly stored, see expiration date for the stability of the kit.

**Safety Information**

Proteinase K Danger

H334



2/6

---

 P304+P341, P342+P311

DTT



Warning

H315, H319

P280, P305+P351+P338, P321, P362, P332+P313,

P337+P313

### Hazard Statements

H315 Causes skin irritation

H319 Causes serious eye irritation

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled

### Precautionary Statements

P280 Wear protective gloves/protective clothing/eye protection/ face protection

P342+P311 If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

P305+P351+ P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

P304+P341 IF INHALED: if breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing

P321 Specific treatment (see on this label)

P362 Take off contaminated clothing and wash before reuse

P332+P313 If skin irritation occurs: Get medical advice/attention

P337+P313 If eye irritation persists: Get medical advice/attention

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information please consult the appropriate material safety data sheets (MSDS).

### Before starting

- **Equilibrate nexttec™ cleanPlates96**

|    |  |
|----|--|
| E1 | Add <b>350 µl Prep Solution</b> onto each well of a nexttec™ cleanPlate96. Incubate for at least <b>5 min</b> at room temperature and centrifuge at <b>350 x g</b> for <b>1 min</b> or apply vacuum for <b>30 to 60 sec</b> to remove excess buffer. |
|----|--|

|    |  |
|----|--|
| E2 | Discard the Waste deep-well plate.<br>Place the <b>nexttec™ cleanPlate96</b> onto a new deep-well plate.<br>Use equilibrated <b>nexttec™ cleanPlates96</b> or store <b>closed at +2 °C to +8 °C</b> and use within <b>one week</b> . |
|----|--|

- Preheat an incubator to **56 °C**

## Protocol

### Lysis

|    |   |
|----|---|
| L1 | Transfer <b>600 µl cold Solution B</b> to each well of a deep-well plate. Add <b>200 µl cold EDTA</b> or <b>heparinized blood</b> . Mix thoroughly by pipetting up and down.  |
| L2 | Incubate deep-well plate <b>5 min on ice</b> . Centrifuge ( <b>2,000 x g, 10 min</b> ), remove <b>700 µl</b> and discard the supernatant.   |
| L3 | Add <b>500 µl cold Solution B</b> to each well of the deep-well plate. Mix thoroughly by pipetting up and down.   |
| L4 | Incubate deep-well plate <b>5 min on ice</b> . Centrifuge ( <b>2,000 x g, 10 min</b> ), remove <b>550 µl</b> and discard the supernatant.   |
| L5 | Add <b>125 µl Buffer R</b> , <b>20 µl Proteinase K</b> and <b>1.5 µl DTT*</b> to each sample. Resuspend pellet by pipetting up and down ( <b>3 times</b> ). Close the deep-well plate using an alu sealing tape. Incubate with shaking ( <b>56 °C, 600 rpm, 30 min</b> ). |

\*For Pre-Mixes see Technical Section.

### Purification of DNA

|   |   |
|---|---|
| P   | Transfer <b>100 µl</b> of the lysate to an <b>equilibrated nexttec™ cleanPlate96</b> .<br>Incubate for <b>3 min at room temperature</b> .<br>Centrifuge at <b>700 x g</b> for <b>1 min</b> or apply vacuum for <b>1 min</b> . |
| <b>The eluate contains the purified DNA!!</b> |   |

### Notes:

**Technical Section****Preparation of Lysis Pre-Mix**

|   |                         |          |         |          |          |
|---|-------------------------|----------|---------|----------|----------|
| LR  | <b>Lysis Buffer LR:</b> | 1 sample | 1 plate | 2 plates | 4 plates |
|   | <b>Buffer R</b>         | 125 µl   | 13.8 ml | 27.5 ml  | 55 ml    |
|   | <b>Proteinase K</b>     | 20 µl    | 2.2 ml  | 4.4 ml   | 8.8 ml   |
|   | <b>DTT</b>              | 1.5 µl   | 165 µl  | 330 µl   | 660 µl   |
| Mix by vortexing. Add <b>146.5 µl</b> of <b>Buffer LR</b> to each sample (L5). The <b>Lysis Buffer LR</b> is stable for <b>1 working day</b> if stored at <b>+2 °C to +8 °C</b> . |                         |          |         |          |          |

- Determination of DNA concentration in nexttec™ 1-Step DNA preparations**

We recommend to determine the DNA concentration:

- Using the fluorescent dye Picogreen® or similar.
- Comparing the fluorescence intensity of DNA bands of unknown concentrations with standards, e.g. in ethidium bromide stained agarose gels.

Please notice:

The use of absorption measurement at 260 nm ( $A_{260}$ ) in a spectrophotometer (e.g. NanoDrop®) for determination of DNA concentration is system related not recommended.

For details and possible workarounds for your specific application please contact:

**service@nexttec.biz**.

- Centrifugation**

For centrifugation at 350x g or 700x g use the settings for relative centrifugal force (RCF) of your centrifuge. Alternatively measure the distance of the nexttec™ cleanPlate96 to the centre of your rotor and calculate the necessary rotations per minute (e.g. rpm =  $299.07 \times \sqrt{350 / r}$  ; r=radius in cm).

- Vacuum Application**

Assemble the vacuum manifold as described in the manual. The flow rate should be between 15 L/min and 150 L/min.

**A regulation of vacuum is not necessary.**

**Product Use Restriction**

nexttec™ 1-Step DNA Isolation Kit components were developed, designed and sold **for research purposes only**. They are suitable for in vitro uses only. No claim or representation is intended for use to identify any specific organism or for clinical use.

It is the responsibility of the user to verify the use of the nexttec™ 1-Step DNA Isolation Kit for a specific application as the performance characteristic of this kit has not been verified to a specific organism.

### Troubleshooting, FAQ and Special Applications

Product claims are subject to change. Therefore, please, visit our website or contact our technical service team for troubleshooting guide, up-to-date protocols and latest applications on nexttec™ 1-Step products.

### Contact Information

Web: [www.nexttec.biz](http://www.nexttec.biz)

Email: [service@nexttec.biz](mailto:service@nexttec.biz)

Tel: +49 (0) 8250 / 92 790 32

Fax: +49 (0) 8250 / 92 790 99

### Ordering Information

For ordering information please visit our website [www.nexttec.biz](http://www.nexttec.biz).

### Distributor

**Biozym**  
SCIENCE IS OUR BUSINESS  
[www.biozym.com](http://www.biozym.com)

#### Biozym Scientific GmbH

Tel.: 05152 / 9020, Fax: 05152 / 2070  
Mail: [support@biozym.com](mailto:support@biozym.com)



#### Biozym Biotech Trading GmbH

Tel.: 01 / 334 0156 0, Fax: 01 / 334 0156 88  
Mail: [support@biozym.com](mailto:support@biozym.com)



6/6