

## User's Instructions

Catalogue number	Size	Scan to product page	Scan to chemical page
62 5121	Extra Small (XS/6)		
62 5122	Small (S/7)		
62 5123	Medium (M/8)		
62 5124	Large (L/9)		
62 5125	Extra Large (XL/10)		

### Determination of material resistance to permeation by chemicals (EN 16523-1:2015) and degradation (EN 374-4:2013)

 ISO 374-1 Type B JKPT	Code letter and chemical	CAS number	Permeation performance and breakthrough time	Degradation (DR) Mean in %
	J: n-Heptane 99%	142-82-5	Level 2, > 30 min.	47.8
	K: Sodium hydroxide 40%	1310-73-2	Level 6, > 480 min.	10.0
	P: Hydrogen peroxide 30%	7722-84-1	Level 6, > 480 min.	35.0
	T: Formaldehyde 37%	50-00-0	Level 6, > 480 min.	25.0

### Marking of gloves protecting against viruses, bacteria and fungi, as defined in ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms – Part 5: Terminology and performance requirements for micro-organisms risks). Protection against viruses: tested according to ISO 16604:2004 (Procedure B) Phi-X, Herpes, HIV.

 ISO 374-5:2016 Level 3 VIRUS	Performance level	Acceptable quality level unit	Inspection level
	Level 3	< 0.65	G1
	Level 2	< 1.5	G1
	Level 1	< 4.0	S4

**(en)** Storage recommendations: Store in a cool and dry place. Gloves should be protected from exposure to direct sunlight, intense artificial light, X-ray machines and other sources of ozone. Disposal: Dispose of in accordance with local authority regulations. This is an MDD Class 1 examination glove (according to Directive 93/42/EEC) and a PPE Category III protective glove according to Regulation (EU) 2016/425. Gloves provide protection against dangerous chemicals and micro-organisms risks. Caution: These gloves should not be used for mechanical and/or thermal protection. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Declaration of conformity accessible from: <http://shieldscientific.com/resource-center/glove-declaration-of-conformity>.

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Website: [www.shieldscientific.com](http://www.shieldscientific.com) ° E-mail: [info@shieldscientific.com](mailto:info@shieldscientific.com) ° ™: trademark of SHIELD Scientific.
All rights reserved. © Copyright 2019.
**EC REP SHIELD Scientific B.V.**
SGS Fimko Oy (Notified Body No: 0598) Takomotie 8, FI-00380 Helsinki, Finland
Biozym Scientific GmbH | Steinbrinksweg 27 | D-31840 Hessisch Oldendorf
SHIELD Scientific B.V. ° Dr. Willem Dreeslaan 1 ° 6721 ND Bennekom ° The Netherlands ° Tel: +31(0)317700202