User's Instructions

Catalogue number	Size	Scan to product page	Scan to chemical page
65 4121	Extra Small (XS/6)		
65 4122	Small (S/7)		
65 4123	Medium (M/8)		
65 4124	Large (L/9)		
65 4125	Extra Large (XL/10)		

Determination of material resistance to permeation by chemicals (EN 16523-1:2015) and degradation (EN 374-4:2013)

 ISO 374-1 Type B KPT	Code letter and chemical	CAS number	Permeation performance and breakthrough time	Degradation (DR) Mean in %
	K: Sodium hydroxide 40%	1310-73-2	Level 6, > 480 min.	-13.3
	P: Hydrogen peroxide 30%	7722-84-1	Level 5, > 240 min.	6.0
	T: Formaldehyde 37%	50-00-0	Level 6, > 480 min.	-15.2

Marking of gloves protecting against viruses, bacteria and fungi, as defined in ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms – Part 5: Terminology and performance requirements for micro-organisms risks). Protection against viruses: tested according to ISO 16604:2004 (Procedure B) Phi-X, Herpes, HIV.

 ISO 374-5:2016 Level 2 VIRUS	Performance level	Acceptable quality level unit	Inspection level
	Level 3	< 0.65	G1
	Level 2	< 1.5	G1
	Level 1	< 4.0	S4

(en) This product contains natural rubber latex
 (fr) Ce produit contient du latex naturel
 (de) Dieses Produkt enthält Naturlatex

(en) This glove contains 50 micrograms or less of total water extractable protein per gram. Caution: This product contains natural rubber latex, which may cause allergic reactions. Warning: This product should not be worn by, or exposed to, individuals allergic to natural rubber latex. Allergy symptoms: If symptoms such as hives, itching, runny nose, watery eyes, facial swelling, abdominal cramps, diarrhoea, nausea, breathing difficulty, rapid heart rate, a sudden drop in blood pressure, dizziness or anaphylactic shock occur, discontinue use of this product. Storage recommendations: Store in a cool and dry place. Gloves should be protected from exposure to direct sunlight, intense artificial light, X-ray machines and other sources of ozone. Disposal: Dispose of in accordance with local authority regulations. This is an MDD Class 1 examination glove (according to Directive 93/42/EEC) and a PPE Category III protective glove according to Regulation (EU) 2016/425. Gloves provide protection against dangerous chemicals and micro-organisms risks. Caution: These gloves should not be used for mechanical and/or thermal protection. These gloves are intended for use in light-duty manufacturing and industrial applications where the demand for the advantages of a shorter glove outweighs the need for additional length. Size 9 (L) and 10 (XL) gloves are shorter in length than required by EN 420:2003 + A1:2009. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Declaration of conformity accessible from: <http://shieldscientific.com/resource-center/glove-declaration-of-conformity>.

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Website: www.shieldscientific.com ° E-mail: info@shieldscientific.com ° ™: trademark of SHIELD Scientific.
All rights reserved. © Copyright 2019.
EC REP SHIELD Scientific B.V.
SGS Fimko Oy (Notified Body No: 0598) Takomotie 8, FI-00380 Helsinki, Finland
Biozym Scientific GmbH | Steinbrinksweg 27 | D-31840 Hessisch Oldendorf
SHIELD Scientific B.V. ° Dr. Willem Dreeslaan 1 ° 6721 ND Bennekom ° The Netherlands ° Tel: +31(0)317700202