



Protocol DNA Purification
Plasmid from Bacteria
with 4ward-NA® puriSorb®

- 4ward-NA® puriPlate96 -

Cat. No. 4w05.P01

Cat. No. 4w05.P04

Cat. No. 4w05.P10

Version 1.3

For research use only

Principle

4ward-NA® DNA purification kits (Plasmid from Bacteria) are designed for extraction and purification of Plasmid DNA from bacteria (high yield/ low yield). **4ward-NA® puriSorb® technology** allows the purification in one single step after lysis. Impurities are bound and the purified DNA is recovered in the eluate. The excellent purity of the DNA allows immediate use in all common downstream applications (e.g. PCR).

Kit contents

The kit contains all necessary reagents for lysis and DNA purification.

Component	Art. No. 4w05.P01	Art. No. 4w05.P04	Art. No. 4w05.P10*
Activation Solution (ActivS)	40 mL	155 mL	400 mL
Buffer 4wE (P)	15 mL	58 mL	144 mL
Protease	0.6 mL	2.25 mL	5.6 mL
RNase A	1.7 mL	6.7 mL	17 mL
Lysozyme (25 mg/mL)**	10.0 mg**	35.0 mg**	85.0 mg**
4ward-NA® puriPlate96	1	4	10
Deep well plate (1.2 mL)	3	12	30
Sealing tapes (plastic)	3	12	30
Storage sealing tapes (alu)	2	8	20

*other kit sizes on request

** add H₂O

Storage Conditions

During shipment all kit components are stable at room temperature. After arrival, store chemicals at **+2 °C to +8 °C**. All other components can be stored at room temperature. Please see expiration date of the kit.

Safety Information

Protease	Danger	H315, H318, H334, H335 P261, P280, P342+P311, P305+P351+P338
----------	--------	-----------------------------------------------------------------

Hazard Statements

H315	Causes skin irritation.
H318	Causes serious eye damage.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

Precautionary Statements

P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P342+P311	If experiencing respiratory symptoms: Call a POISON CENTER or doctor/ physician.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please contact us for the appropriate material safety data sheets (MSDS).

Before starting

Activation of puriSorb® in 4ward-NA® puriPlate96

A1	Add 350 µL Activation Solution onto each well of a 4ward-NA® puriPlate96 . Incubate for at least 5 min at room temperature. Centrifuge at <u>max.</u> 350 × g for 1 min or apply vacuum for 1 min to remove excess Activation Solution. Discard the deep well plate.
A2	The activated 4ward-NA® puriPlate96 is ready for purification.

- **Preheat an incubator to 96 °C.**
- **Dissolve Lysozyme**
Add **0.4 mL** (4w03.P01), **1.4 mL** (4w03.P04) or **3.4 mL** (4w03.P10) of sterile water to the lyophilized lysozyme (final concentration is 25 mg/mL).

Protocol

Lysis of Bacteria (Plasmid) for DNA purification

1. Transfer **up to 1 mL of bacterial culture** into each well of a deep well plate (1.2 mL).
Centrifuge (**2,000 × g, 10 min**) and remove carefully the supernatant.
2. Add **130 µL Buffer 4wE (P)**, **5 µL Protease** and **3 µL Lysozyme** to the cell pellet.
Dissolve cell pellets and close plate with an alu sealing tape.
Incubate for **5 min at room temperature**.
3. Transfer plate to a preheated incubator (**96 °C**) and incubate for **10 min**.
Let the plate cool to room temperature.
4. Add **15 µL RNase A** to each lysate. Mix gently. Incubate for **3 min** at room temperature.

Purification of DNA using 4ward-NA® puriSorb® (ONE-STEP)

P1	Transfer between 60 µL and 120 µL of the lysate into an activated well of the 4ward-NA® puriPlate96 . Incubate for at least 3 min and centrifuge at <u>max.</u> 700 × g for 1 min .
	The eluate contains the ready-to-use purified DNA.

Troubleshooting, FAQ and Special Applications

Please contact us for further information about our innovative DNA Purification System using the **4ward-NA® puriSorb®** technology.

Contact Information

Web: www.4ward-na.com
Email: service@4ward-na.com
Tel: +49 175 948 33 46

Ordering Information and Prices

For ordering information and prices please contact: sales@4ward-na.com.

