

# Biozym HS Taq DNA Polymerase

## Product information

Store at -20°C.

The kit may be stored at 4°C for short term use (up to 1 month). Multiple freeze/thaw cycles (up to 30) are possible with no loss of activity.

Component	331620X 50 Units (5 U/μl)	331620L 500 Units (5 U/μl)	331620XL 2500 Units (5 U/μl)
Biozym HS Taq DNA Polymerase	1 x 5 μl	2 x 50 μl	10 x 50 μl
5 x Reaction Buffer	1 x 0.4 ml	4 x 1 ml	20 x 1 ml

For research use only.

## 1. Introduction

Biozym HS Taq DNA Polymerase is designed for use in all routine and high throughput PCR applications including multiplex PCR, genotyping, screening, library construction, colony PCR and direct PCR. The robust enzyme ensures consistent results on a broad range of templates.

Using latest technology the polymerase enables short extension times, enhanced specificity and improved yields – ideal for complex templates and multiplex PCR. The used hot-start technology inhibits DNA polymerase activity at ambient temperature, thus preventing amplification of non-specific products.

## 2. Notes

- Biozym HS Taq DNA Polymerase produces A-tailed ends.
- Error rate: Approximately 1 error per  $2.0 \times 10^5$  nucleotides incorporated.
- Biozym HS Taq DNA Polymerase is particularly resistant to PCR inhibitors. Direct PCR could be performed from unprocessed samples including bacterial culture, bacterial colonies, blood and urine. (see 5.8)

## 3. Technical support

For technical support please contact [support@biozym.com](mailto:support@biozym.com)

### Biozym Scientific GmbH

Steinbrinksweg 27  
D-31840 Hess. Oldendorf  
Ph: +49 5152 9020  
[www.biozym.com](http://www.biozym.com)

#### 4. Reaction conditions for PCR

Carefully mix and centrifuge the tube before opening to ensure homogeneity and improve recovery.

**Table 1. Pipetting instructions**

Component	50 µl reaction	20 µl reaction	Final conc.	Notes
5x Reaction Buffer	10 µl	4 µl	1x	
Biozym HS Taq DNA Polymerase	0.25 – 1.0 µl	0.1 – 0.4 µ		
Forward primer (10 µM)	2.0 µl	0.8 µl	400 nM	See 5.2 for primer design
Reverse primer (10 µM)	2.0 µl	0.8 µl	400 nM	
Template DNA cDNA gDNA	<100 ng <0.5 µg	<40 ng <0.2 µg	Variable	See 5.3 for template
PCR grade water	Up to 50 µl reaction volume	Up to 20 µl reaction volume		

**Table 2. Cycling instructions**

Cycles	Temperature	Time	Notes
1	95°C	1 - 2 min	Initial denaturation and enzyme activation. Colony PCR: 10 min.
40	95°C 55 to 65°C 72 °C	15 seconds 15 seconds 1 - 90 seconds	Denaturation Annealing Extension (15 s/kb). For multiplex PCR use 90 seconds

#### 5. Notes about reaction components and cycling conditions

##### 5.1. Reaction Buffer

The 5x reaction buffer contains: 15 mM MgCl<sub>2</sub>, 5 mM dNTPs, enhancers and stabilizers. We do not recommend adding further PCR enhancers or MgCl<sub>2</sub>.

##### 5.2. Primers

Primers should have a predicted melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>). The final primer concentration in the reaction should be between 0.2 µM and 0.6 µM.

##### 5.3. Template

For eukaryotic DNA use between 5 ng and 500 ng per reaction, for cDNA use below 100 ng per reaction.

##### 5.4. Annealing

Use an annealing temperature equal to the T<sub>m</sub> of the lower T<sub>m</sub> primer. If non-specific products are present increase in 2°C increments. Alternatively use a temperature gradient to experimentally find the optimal annealing temperature.

##### 5.5. Extension

Extension should be performed at 72°C. The optimal extension time is dependent on amplicon length and complexity of template. We recommend an extension time of 15 seconds per kilobase (kb) for amplification from eukaryotic DNA (for amplicons between 1 kb and 6 kb).

##### 5.6. Multiplex PCR

When first performing multiplex PCR it is recommended to run a temperature gradient around the calculated annealing temperatures. The annealing temperature representing best specificity should be used in subsequent experiments. Fast cycling conditions should not be used. Initially a 90 second extension time is recommended. To increase yield this time may be extended further.

##### 5.7. Colony PCR

Use a sterile tip to pick a bacterial colony and resuspend into a 50µl reaction mix as described in table 1 above. From liquid culture add 5µl of overnight culture to the master mix. For lysis of cells initial denaturation time should be 10 minutes.

##### 5.8. Direct blood/urine PCR

Add 2µl mammalian blood or urine to a 50µl reaction mix as described in table 1.