

## Package insert

### S7 Fusion Polymerase™

V2-0

Product no. MD-S7-100 / MD-S7-500

Biozym Art. nr.: 332530S, 332530L

#### 1. Introduction

S7 Fusion Polymerase™ is a high-fidelity DNA polymerase that offers extraordinary performance in PCR applications. The polymerase is fused to the processivity enhancing Sso7d domain from *Sulfolobus solfataricus*, making S7 Fusion Polymerase™ one of the fastest and most robust PCR polymerases on the market. The polymerase possesses a 3' → 5' exonuclease activity and it generates blunt ends in PCR.

#### 2. Reagent kit components

Component name	Product number	MD-S7-100	MD-S7-500
S7 Fusion Polymerase™ 2 U/μl	MD-E604S	100 U, 50 μl	-
S7 Fusion Polymerase™ 2 U/μl	MD-E604L	-	500 U, 250 μl
5X HF buffer	MD-B704	2 x 1,5 ml	6 x 1,5 ml
5X GC buffer	MD-B703	1,5 ml	2 x 1,5 ml
50 mM MgCl <sub>2</sub>	MD-M501	1,5 ml	2 x 1,5 ml
DMSO	MD-D502	500 μl	500 μl

#### 3. Shipping and storage

The S7 Fusion Polymerase™ reagents are packed into a plastic pouch and shipped on gel ice. After reception, the product should be stored at -20°C.

If stored under the recommended conditions the product will maintain its performance until the expiry date.

#### 4. Guidelines for using S7 Fusion Polymerase™

##### PCR setup

Prepare the PCR reactions according to the table below. Add S7 Fusion Polymerase™ last.

	20 μl reaction	50 μl reaction	Final concentration
Water	Add to 20 μl	Add to 50 μl	
5X HF or 5X GC buffer	4 μl	10 μl	1X
10 mM dNTPs	0,4 μl	1 μl	200 μM each
Forward primer	X μl	X μl	0,5 μM
Reverse primer	X μl	X μl	0,5 μM
Template DNA	See instruction for details		
S7 Fusion Polymerase™	0,2 μl	0,5 μl	0,02 U/μl

##### Buffers and additives

Two buffers are provided with the S7 Fusion Polymerase™. The HF buffer is recommended for high-fidelity applications, as the error rate in PCR is lower than with the GC buffer. The GC buffer is recommended for GC-rich templates or long PCR reactions. The GC buffer may contain some precipitate after thawing; the precipitate dissolves upon heating.

GC-rich or long PCR reactions might need DMSO in order to produce a satisfactory yield. Add 3% DMSO (e.g. 1,5 μl per 50 μl reaction) as a starting point for DMSO optimization. High concentrations of DMSO affect the T<sub>m</sub> of the primers. As a general guideline, a final concentration of 10% DMSO requires a 5°C decrease in the annealing temperature. Other PCR additives such as formamide, glycerol and betaine can also be used for difficult templates.

##### MgCl<sub>2</sub>

The 5X HF and 5X GC buffers contain an optimal concentration of MgCl<sub>2</sub>. Higher Mg<sup>2+</sup> concentrations can increase PCR yield but decrease specificity, whereas low Mg<sup>2+</sup> concentrations may lead to low or no yield. An increase of the MgCl<sub>2</sub> concentration may be needed if the template or primers contain Mg<sup>2+</sup> binding chelators, such as EDTA or EGTA. Also, if there is an excess of dNTPs present in the reaction (e.g. when using cDNA synthesis mix as template), the amount of MgCl<sub>2</sub> should be increased. When adjusting the MgCl<sub>2</sub> concentration, make the increments in 0,2 mM steps (e.g. 0,2 μl 50 mM MgCl<sub>2</sub> per 50 μl reaction).

## dNTPs

It is recommended to use 200 µM of each dNTP. NOTE: Do not use dUTP, since S7 Fusion Polymerase™ cannot amplify templates containing dUTP.

## Template DNA

The amount of needed template DNA is dependent on the complexity. For low complexity templates, such as plasmid or lambda DNA, a few picograms to a few nanograms of template is sufficient in a 50 µl PCR reaction. For high complexity genomic templates, 50-250 ng per 50 µl PCR reaction can be used. If using a cDNA synthesis reaction mixture as template, increasing the template amount to more than 10% of the PCR mixture does generally not increase yield.

## Polymerase

It is recommended to use 1 unit S7 Fusion Polymerase™ per 50 µl PCR reaction, but the optimum can vary between 0,5 and 2 units per 50 µl PCR reaction. For amplicons that are several kb long, no more than 2 units per 50 µl should be used.

## PCR cycling

For optimal performance, the following PCR cycling protocols are recommended:

### 2-step protocol

Cycles	Step	Temperature	Time
1	Initial denaturation	98°C	30 s
25-35	Denaturation	98°C	5-10 s
	Annealing and extension	72°C	15-30 s/kb
1	Final extension	72°C	5-10 min
1	Cooling	4°C	hold

### 3-step protocol

Cycles	Step	Temperature	Time
1	Initial denaturation	98°C	30 s
25-35	Denaturation	98°C	5-10 s
	Annealing	X°C	10-30 s
	Extension	72°C	15-30 s/kb
1	Final extension	72°C	5-10 min
1	Cooling	4°C	hold

## Denaturation

A short denaturation time of 5-10 seconds at 98°C is sufficient for most templates.

## Annealing

The S7 Fusion Polymerase™ requires higher annealing temperatures than most other PCR polymerases. The higher annealing temperature allows for more frequent use of a two-step protocol, where the annealing and extension steps are combined. As a rule of thumb, for primers >20 nt in length, the annealing temperature should be the T<sub>m</sub> of the lower T<sub>m</sub> primer +3°C. For primers ≤20 nt, the annealing temperature should be the T<sub>m</sub> of the lower T<sub>m</sub> primer. The two-step protocol should be used when the annealing temperatures of the primers are 72°C or higher.

## Extension

Amplicon complexity and length determine the required extension time. Simple templates require 15 s/kb extension time, whereas more complex templates need 30 s/kb.

## 5. Specifications

### Storage buffer

20 mM Tris-HCl pH 8,0, 100 mM KCl, 0,1 mM EDTA, 0,2 mg/ml BSA, stabilizers and 50% glycerol.

### Unit definition

One unit is defined as the amount of enzyme that will incorporate 10 nmoles of dNTPs into a polynucleotide fraction at 72°C in 30 min.

### Endonuclease contamination assay

No conversion to nicked form is detected after incubating S7 Fusion Polymerase™ with supercoiled ΦX174 DNA.

### Exonuclease contamination assay

No contaminating exonuclease activity is detected after incubating S7 Fusion Polymerase™ with radiolabeled DNA.

### DNA amplification

S7 Fusion Polymerase™ amplifies a 7,5 kb genomic DNA fragment and a 10 kb lambda DNA fragment.

### Safety

This product should be handled only by persons trained in laboratory techniques.

### Terms and conditions of sale

All goods and services are sold subject to the terms and conditions of sale of duly authorized officer of Mobidiag Ltd. A copy of these terms and conditions is available upon request. Prices are subject to change without prior notice.

### Notice to the user

The information presented here is accurate and reliable in the best of our knowledge and belief, but is not guaranteed to be so. The information in this document is of a proprietary nature. Mobidiag is the copyright holder of the Mobidiag materials. No part of this document may be reproduced, copied, or transmitted in any form or by any means without a prior written permission of Mobidiag. Nothing herein is to be construed as recommending any practice or any product in violation of any patent or in violation of any law or regulation.

This product is licensed from Bio-Rad Laboratories, Inc. under U.S. Pat. Nos. 6,627,424, 7,541,170, 7,670,808, 7,919,296, 8,232,078, 8,415,129, 8,476,045, 8,895,283, 8,900,846, and 9,139,873, and 7,560,260, 8,367,376, and 8,470,573, and related patents including corresponding patents in other countries, for use only in genomic applications in standard (non-real-time) PCR in the life science research field. The purchaser/end-user of this product is not licensed or authorized to use the product in any in-vitro diagnostic, industrial microbiology, real-time PCR, or digital PCR application.

Mobidiag is an ISO 9001 and ISO 13485 certified company.

© Copyright, Mobidiag. All rights reserved.  
Document version 2-1