CapClip Pyrophosphatase 5 U/ μ l

Cellscript

CELLSCRIPT™
RNA for Translation in Cells

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Cap-Clip™ Acid Pyrophosphatase (100 units)

20 μ l Cap-Clip™ Acid Pyrophosphatase (5 U/ μ l) in 50% Glycerol, 10 mM Tris-HCl, pH 7.5, 100 mM NaCl, 1 mM Dithiothreitol, 0.1 mM EDTA and 0.01% Triton® X-100.

250 μ l 10X Cap-Clip™ Acid Pyrophosphatase Reaction Buffer (0.5 M NaOAc, pH 6.0, 10 mM EDTA, 1% Beta-Mercaptoethanol and 0.1% Triton X-100)

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Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

CELLSCRIPT and Cap-Clip are trademarks of CELLSRIPT, Madison, Wisconsin.

© 2015 CELLSRIPT, All rights reserved.

726 Post

Cap-Clip™ Acid Pyrophosphatase is a plant-derived enzyme that hydrolyzes various pyrophosphate bonds, including the pyrophosphate bonds of the 5'-terminal m7GpppG "cap" of eukaryotic messenger RNAs, as well as 5' cap structures on many small nuclear RNAs (snRNAs), heterogeneous nuclear RNAs and some viral RNAs. Complete hydrolysis of such capped RNAs generates RNA that has a 5'-monophosphate group.

Applications

- Cloning of Full-Length cDNAs and Rapid Amplification of cDNA Ends (RACE): After decapping the RNA an adaptor sequence is ligated to the 5' end using T4 RNA
- RNA 5'-End Conversion. In addition to decapping capped mRNA, TAP can also be used to convert 5'-triphosphate RNA into 5'-monophosphate RNA.

MATERIALS

Materials Supplied

Store at -20°C in a freezer without a defrost cycle. Do not store at -70°C.

Cap-Clip™ Acid Pyrophosphatase (100 units)

Cap-Clip™ Acid Pyrophosphatase, 5 U/μl in 50% glycerol, 10 mM Tris-HCl, pH 7.5, 100 mM NaCl, 1 mM dithiothreitol, 0.1 mM EDTA and 0.01% Triton® X-100.

20 μl

10X Cap-Clip™ Acid Pyrophosphatase Reaction Buffer 0.5 M NaOAc, pH 6.0, 10 mM EDTA, 1% beta-mercaptoethanol and 0.1% Triton X-100.

250 μl

SPECIFICATIONS

Unit Definition

One unit of Cap-Clip Acid Pyrophosphatase releases one nanomole of inorganic phosphate from m7GpppG in 30 minutes at 37°C under standard assay conditions.

Contaminating Activity Assays

Cap-Clip Acid Pyrophosphatase is free of detectable RNase and DNase activities.

Art.Nr Art.Bez.1 Art.Bez.2 Preis

187005 CapClip Pyrophosphatase 5 U/μl 100 Units 475,00 €

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

CELLSCRIPT and Cap-Clip are trademarks of CELLSRIPT, Madison, Wisconsin.

© 2015 CELLSRIPT, All rights reserved.

726 Post