

Kammoptionen für die horizontale Gelelektrophorese

EasyPhor Mini				EasyPhor Midi				EasyPhor Medi				EasyPhor Maxi			
Artikel Nr.	Anzahl Zähne	Stärke (mm)	Volumen (µl)	Artikel Nr.	Anzahl Zähne	Stärke (mm)	Volumen (µl)	Artikel Nr.	Anzahl Zähne	Stärke (mm)	Volumen (µl)	Artikel Nr.	Anzahl Zähne	Stärke (mm)	Volumen (µl)
615021	1+1	0.75	152	615099	1+1	0.75	270	615182	1+1	0.75	371	615272	1+1	0.75	506
615022	2+2	0.75	68	615100	2+2	0.75	118	615183	2+2	0.75	169	615273	2+2	0.75	236
615023	4+2	0.75	36	615101	4+2	0.75	57	615184	4+2	0.75	91	615274	4+2	0.75	115
615024	8 MC	0.75	8	615102	8	0.75	30	615185	10	0.75	34	615275	10	0.75	54
615025	8	0.75	19	615103	10 MC	0.75	20	615186	12	0.75	30	615276	16	0.75	30
615026	10	0.75	14	615104	12	0.75	17	615187	20	0.75	16	615277	20 MC	0.75	20
615027	12 MC	0.75	10	615105	16	0.75	12	615188	35	0.75	7	615278	25	0.75	16
615028	16	0.75	7	615106	20 MC	0.75	10	615189	10 MC	0.75	22	615279	30	0.75	13
615029	1+1	1	203	615107	25	0.75	7	615190	14 MC	0.75	22	615280	36	0.75	11
615030	2+2	1	90	615108	1+1	1	360	615191	16 MC	0.75	20	615281	40 MC	0.75	8
615031	4+2	1	48	615109	2+2	1	158	615192	18 MC	0.75	8	615282	50	0.75	8
615032	8 MC	1	11	615110	4+2	1	77	615193	28 MC	0.75	8	615283	1+1	1	675
615033	8	1	25	615111	8	1	41	615194	30 MC	0.75	9	615284	2+2	1	315
615034	10	1	18	615112	10 MC	1	27	615195	1+1	1	495	615285	4+2	1	153
615035	12 MC	1	14	615113	12	1	23	615196	2+2	1	225	615286	10	1	72
615036	16	1	10	615114	16	1	16	615197	4+2	1	122	615287	16	1	41
615037	1+1	1.5	304	615115	20 MC	1	14	615198	10	1	45	615288	20 MC	1	27
615038	2+2	1.5	135	615116	25	1	10	615199	12	1	41	615289	25	1	21
615039	4+2	1.5	72	615117	1+1	1.5	540	615200	20	1	21	615290	30	1	17
615040	8 MC	1.5	17	615118	2+2	1.5	236	615201	35	1	10	615291	36	1	14
615041	8	1.5	37	615119	4+2	1.5	115	615202	10 MC	1	29	615292	40 MC	1	11
615042	10	1.5	27	615120	8	1.5	61	615203	14 MC	1	29	615293	50	1	10
615043	12 MC	1.5	20	615121	10 MC	1.5	41	615204	16 MC	1	27	615294	1+1	1.5	1013
615044	16	1.5	15	615122	12	1.5	34	615205	18 MC	1	11	615295	2+2	1.5	473
615045	1+1	2	405	615123	16	1.5	24	615206	28 MC	1	11	615296	4+2	1.5	230
615046	2+2	2	180	615124	20 MC	1.5	20	615207	30 MC	1	13	615297	10	1.5	108
615047	4+2	2	96	615125	25	1.5	15	615208	1+1	1.5	743	615298	16	1.5	61
615048	8 MC	2	23	615126	1+1	2	720	615209	2+2	1.5	338	615299	20 MC	1.5	41
615049	8	2	50	615127	2+2	2	315	615210	4+2	1.5	182	615300	25	1.5	32
615050	10	2	36	615128	4+2	2	153	615211	10	1.5	68	615301	30	1.5	26
615051	12 MC	2	27	615129	8	2	81	615212	12	1.5	61	615302	36	1.5	22
615052	16	2	20	615130	10 MC	2	54	615213	20	1.5	32	615303	40 MC	1.5	17
				615131	12	2	45	615214	35	1.5	15	615304	50	1.5	16
				615132	16	2	32	615215	30 MC	1.5	19	615305	1+1	2	1350
				615133	20 MC	2	27	615216	10 MC	1.5	44	615306	2+2	2	630
				615134	25	2	20	615217	14 MC	1.5	44	615307	4+2	2	306
								615218	16 MC	1.5	41	615308	10	2	144
								615219	18 MC	1.5	17	615309	16	2	81
								615220	28 MC	1.5	17	615310	20 MC	2	54
								615221	1+1	2	990	615311	25	2	42
								615222	2+2	2	450	615312	30	2	34
								615223	4+2	2	243	615313	36	2	29
								615224	10	2	90	615314	40 MC	2	23
								615225	12	2	81	615315	50	2	21
								615226	16 MC	2	54				
								615227	20	2	43				
								615228	30 MC	2	25				
								615229	10 MC	2	59				
								615230	14 MC	2	59				
								615231	18 MC	2	23				
								615232	28 MC	2	23				
								615233	35	2	20				

Biozym



Eingetragene Warenzeichen: Alpha Innotech Corp.: Alphamager, FluorChem. National Diagnostics Inc.: AcrylaGel, Bis-AcrylaGel

Biozym Scientific GmbH Postfach, D-31833 Hess. Oldendorf, Tel.: 0 51 52 / 90 20, Fax: 0 51 52 / 20 70

Biozym Biotech Trading GmbH Wehlistraße 27 b, A-1200 Wien, Tel.: 01 / 334 0156 0, Fax: 01 / 334 0156 88

www.biozym.com