

# PCR Disposables

## Thermocycler Kompatibilitätstabelle für 96 und 384 Well PCR Platten

Well Formen:		Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.
Standard Well (High Profile), Bauhöhe ca. 21 mm, max. Kapazität > 350 µl		710884	710876	710881	710880	710883	710882	710875	710885
Low Profile Well, Bauhöhe ca. 16 mm, max. Kapazität 200 µl									
Hersteller	Thermocycler Modell	Unskirted High Profile	Unskirted Low Profile	Unskirted Chimney Top High Profile	Semi Skirted High Profile	Semi Skirted Low Profile	Subskirted High Profile	Full Skirted High Profile	Full Skirted 384 Well Platte
Applied Biosystems	PRISM 3100 & 3130(XL)	x		x	x		x		
	PRISM 3700	x		x	x		x		x
	PRISM 3730(XL)	x		x	x		x		
	GeneAmp 2700	x		x	x		x		x
	GeneAmp 9600	x		x	x		x		
	GeneAmp 9700	x		x	x		x		x
	GeneAmp 5700	x			x		x		
	PRISM 7000	x			x		x		
	7300	x			x	x	x		
	7500	x			x	x	x		
7700	x			x		x			
7900HT	x			x	x	x			
Amersham BioScience	MegaBace 500							x	
	MegaBace 1000							x	
	MegaBace 4000								x
Biometra	Uno	x	x	x	x			x	
	Uno II	x	x	x	x				x
	T1	x	x	x	x			x	x
	TGradient	x	x	x	x			x	
	TRobot	x	x	x	x			x	x
BioRad	Gene Cyclyer	x		x	x		x		
	iCycler	x		x				x	
	MyCycler	x		x					
	PTC 100	x	x	x				x	x
	PTC 200	x	x	x				x	x
	PTC 225	x	x	x				x	x
	Dyad und Dyad disciple	x	x	x				x	x
	Opticon/Opticon 2	x	x					x	
	Chromo 4		x					x	
	MyiQ	x		x				x	
Corbett Research	IQ5	x		x				x	
	PalmCycler 96		x					x	
Eppendorf	PalmCycler 384								x
	Mastercycler	x	x	x				x	
	Mastercycler Gradient	x		x	x		x		
	Mastercycler ep	x	x	x				x	
	Mastercycler M384								x
Mastercycler ep Realplex		x					x		
Ericom	Power Block I	x	x	x	x	x	x		
	Deltacycler I	x	x	x					
	Deltacycler II	x	x	x	x	x	x		
	Single Block	x	x	x					
	Twin Block	x	x	x					
G-Storm	GS1	x	x	x	x	x	x		
	GS4	x	x	x	x	x	x		
	GSX	x	x	x	x	x	x		
	GSXs	x	x	x	x	x	x		
MWG	Primus 96	x	x	x	x	x	x	x	
	Primus 384								x
Stratagene	Robocycler 96	x							
	Robocycler Gradient	x	x	x	x	x	x	x	
	MX4000 und MX3000P	x	x	x	x	x	x		
TaKaRa	TP 240							x	
	TP 3000	x	x	x	x	x	x	x	
Techne	Touchgene	x	x	x	x	x	x		
	Cyclogene	x	x	x	x	x	x		
	Genius	x	x	x	x	x	x		
	Genius Quad	x	x	x	x	x	x		
	Genius(TC412)	x	x	x	x	x	x	x	
	Flexigene	x	x	x	x	x	x	x	x
	Touchgene X	x						x	x
	Touchgene Gradient (TC512)	x	x	x	x	x	x	x	x
	Quantica		x					x	
	Hybaid	Omnigene	x	x	x	x	x	x	x
PCR Express & Omni-E		x	x	x	x	x	x	x	x
Px2 & PxE		x	x	x	x	x	x	x	x
MBS Satellite (Multiblock) System		x	x	x	x	x	x	x	x
PCR Sprint		x	x	x	x	x	x	x	x
Touchdown		x	x	x	x	x	x	x	x