



SuperClear[®] Specimen Collection and Transport Tubes

Complete Safety in Specimen Shipment

✓ 95kPa certified as designated by the IATA and U.S. DOT Standard

Find out more at labcon.com

labcon[®]

SuperClear® Specimen Collection and Transport Tubes

Complete Safety In Specimen Shipment

Labcon's SuperClear® Specimen Collection and Transport Tubes provide certified safety when storing and shipping biological material. These leakproof freestanding tubes meet the IATA and U.S. DOT 95kPa performance specifications legally required when shipping samples by ground or air. This ensures transport tubes are leakproof under differential pressure conditions and extreme temperature shifts as outlined in IATA DGR UN3373 and U.S. DOT HMR, 49 C.F.R., Parts 171-180.

Made from 100% virgin USP Class-VI medical-grade resins, SuperClear® Specimen Collection and Transport Tubes are suitable for use in medical, diagnostic and research labs and are certified contamination-free as mandated by the U.S. CDC. These certified sterile tubes have a Sterility Assurance Level (SAL) of 10^{-6} , as required for diagnostic assays.

- ✓ **95kPa certified**
Meets legal requirements for IATA/
U.S. DOT/UN3373 air & ground shipping
- ✓ **Leakproof**
Leakproof under differential pressure
& extreme temperature shifts
- ✓ **Validated Sterility**
 10^{-6} SAL, as required for diagnostic
& medical assays

Product Features

Certified to IATA/U.S. DOT/
UN3373 standards for
domestic and international
ground and air transport

Graduations at 1mL
increments and a
frosted writing surface

Made from non-
cytotoxic medical-
grade virgin resins

Leakproof caps feature
molded-in sealing rings
that form a certified
barrier against air, liquid,
and aerosol leaks

Designed for secure
shipment of viral
transport media and
biological samples

Freestanding base
includes a concentrator
cone for centrifuging
and holding swabs

Quality & Compliance



5 mL SuperClear® Specimen Collection and Transport Tubes

Part Number	Sterile	Packaging	Pk	Case
LC0400N		Bulk, Tubes & Caps	500	1.000
LC0403N	✓	Bulk, only Tubes	500	1.000
LC0404N	✓	Bulk, only Caps	500	1.000

10 mL SuperClear® Specimen Collection and Transport Tubes

Part Number	Sterile	Packaging	Pk	Case
LC0410N		Bulk, Tubes & Caps	500	1.000
LC0413N	✓	Bulk, only Tubes	500	1.000
LC0404N	✓	Bulk, only Caps	500	1.000

Contact us today:

Biozym
SCIENCE IS OUR BUSINESS
www.biozym.com

Biozym Scientific GmbH
Tel.: +49 (0) 5152 9020
Mail: support@biozym.com

Biozym Biotech Trading GmbH
Tel.: +43 (0) 1 334 0156 0
Mail: support@biozym.com

