



### Squid™ Pipette

Full Range Pipette (1 - 1000 µL)

#### Simplify Your Bench

- Replaces a Full Set of Pipettes
- 1 - 1000 µL in a Single Device
- Universal Tip Compatibility
- One-Click Switching Between Tip Sizes

#### How Does It Work?

The unrivaled performance of Squid is achieved through **Dynamic Volume Control™**, which combines a patented tip selector mechanism and DeNovix expertise in microvolume motor control.

Squid employs two tip sizes, spanning the full 1- 1000 µL range to eliminate the need to switch devices mid-workflow. Instantly activate each mode with the single press of a button.

Switch with the Press of a Button



**Large Tip Mode**  
100 - 1000 µL

**Small Tip Mode**  
1 - 200 µL

#### Designed by Life Scientists

- **Ergonomic Design:** A comfortable, balanced grip for prolonged use
- **Custom Presets:** Instantly switch between required volumes
- **Lab-Ready:** Pre-configured modes and custom options available
- **All-Day Battery Life:** Equipped to support your workflow

#### Exceeds ISO Calibration Specifications Across the Full Range

The Squid Pipette outperforms industry standards across its entire 1 - 1000 µL range, doing the work of four traditional pipettes without compromising on data quality.

Squid Performance Against ISO 8655		
ISO Range	Systematic Error (Accuracy)	Random Error (Precision)
1000 - 100 µL	✓	✓
200 - 20 µL	✓	✓
20 - 2 µL	✓	✓
10 - 1 µL	✓	✓



**View Performance Data**

Tested and certified to outperform ISO 8655 requirements from 1 µL to 1000 µL.

## A Full Set in One Pipette

The **DeNovix Squid™** is a revolutionary electronic pipette that replaces an entire traditional set – from **P10 to P1000** – with a single, high-performance device.

# Squid™

Full Range Pipette

Personalize  
Your Lab  
with 5 Color  
Choices



Art. Nr. 315Q-B

Art. Nr. 315Q-R

Art. Nr. 315Q-L

Art. Nr. 315Q-G

Art. Nr. 315Q-D

## Specifications



<b>Pipetting Type</b>	Air Displacement Pipette
<b>Volume Range</b>	1 - 1000 µL
<b>Accuracy/Precision</b>	Exceeds ISO 8655-2_2022 at All Volumes
<b>Warranty</b>	2 Years
<b>Weight</b>	230 grams
<b>Autoclavable</b>	Yes – Tip Holders and Ejectors
<b>Dimensions (LxW)</b>	30.4 x 4.4 cm
<b>Display</b>	1.5 Inch Full-Color HD Screen
<b>Connectivity</b>	Bluetooth 5.0, Wi-Fi 2.4 GHz
<b>Power Supply Input</b>	100 - 240 V ~ 50/60 Hz
<b>Power Supply Output</b>	5 V / 3.0 A / 15.0 W
<b>Battery Type</b>	1200 mAh Li Ion UL-1642 Compliant
<b>Charging</b>	USB-C Cable (Included) Squid Charging Stand (Optional)
<b>Approvals</b>	CE, UL/CSA, FCC and Japan CAB/VCCI
<b>Patents</b>	denovix.com/patents

# Biozym

SCIENCE IS OUR BUSINESS

[www.biozym.com](http://www.biozym.com)

**Biozym Scientific GmbH**

Tel.: +49 (0) 5152 9020  
Mail: [support@biozym.com](mailto:support@biozym.com)



**Biozym Biotech Trading GmbH**

Tel.: +43 (0) 1 334 0156 0  
Mail: [support@biozym.com](mailto:support@biozym.com)