



CE

- + High capacity, benchtop model
- + Circular 19 mm orbit for aeration and mixing
- + Instantly exchange flask clamps with MAGIC Clamp™ accessories

Technical Data:

|                         |                              |
|-------------------------|------------------------------|
| Speed Range:            | Variable, 30 to 300 rpm      |
| Speed Increment:        | 1 rpm                        |
| Motion / Orbit:         | Orbital / 19 mm              |
| Load Capacity:          | 5 x 2 l or 9 x 1 l / (30 kg) |
| Platform Size:          | 45 x 45 cm                   |
| Timer:                  | 1 min. to 48 hr / continuous |
| Operating Temp. Range:  | +4 °C to +65 °C              |
| Dimensions (W x D x H): | 147 x 54 x 12.5 cm           |
| Weight:                 | 29 kg                        |
| Electrical:             | 230V, 50-60Hz                |
| Warranty:               | 2 Years                      |

Ordering Information:

| Item No.     | Description:  |
|--------------|---|
| 55BT1010-E   | Orbi-Shaker XL with flat mat platform (47 x 54 cm)    |
| 55H1010-MR** | MAGic Clamp™ platform for flasks & tubes (47 x 54 cm) |
| 55H1010-P-SP | Universal spring platform for unique vessel types     |
| 55H1010-P-MP | Dedicated platform for 6 x microplates                |

\*\*Flask clamps and tube racks sold separately, see pg. 25 for a complete listing.



CE

- + Digital speed control up to 1,500rpm
- + 3mm orbit for thorough mixing in microplates
- + Includes platform for up to four microplates
- + Mix up to 80 microtubes with optional adapter

Technical Data:

|                         |  |
|-------------------------|--|
| Speed Range:            | 200 to 1,500 rpm                             |
| Orbit:                  | 3 mm   |
| Max. Capacity:          | Up to 4 microplates                          |
| Operating Temp. Range:  | +4 °C to +65 °C                              |
| Dimensions (W x D x H): | 33.3 x 29.8 x 12 cm                          |
| Timer:                  | 1 min. to 99 hour / continuous               |
| Weight:                 | 5.5 kg                                       |
| Electrical:             | Universal Voltage<br>100-240V, 50-60Hz, 250W |
| Warranty:               | 2 Years                                      |

Ordering Information:

| Item No.       | Description:   |
|----------------|--|
| 55BT1502-E     | Orbi-Shaker™ MP microplate shaker, includes platform for 4 microplates |
| 55BT1500-A1520 | Optional adapter, 20 x 1.5/2.0 ml for microplate shakers               |
| 55BT1500-A05   | Optional adapter, 20 x 0.5 ml for microplate shakers                   |