


Chemical resistance guide

LEVEL 0	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	LEVEL 6
< 10 min	10 > 30 min	30 > 60 min	60 > 120 min	120 > 240 min	240 > 480 min	> 480 min

duoSHIELD* PFT Nitrile 290



- Virus resistant / Chemotherapy drugs
- Length: 290 mm/ 11.4"
- Palm thickness: 0.08 mm/ 3.1 mil
- Color: #1073b6
- Design: Ambidextrous / Powder-free
- Allergies: Latex-free / Free of Thiazoles
- Chemical performance: Type B / High chemical protection
- Biological risk:
- Mechanical risk: N/A
- Particles level: N/A / N/A
- Applications: Laboratory / General workplace / Healthcare

50-00-0 Formaldehyde 37%	LEVEL 4 121 min
7722-84-1 Hydrogen peroxide 30%	LEVEL 6 480 min
1310-73-2 Sodium Hydroxide 40%	LEVEL 6 480 min

DISCLAIMER: The data provided was based on gloves tested under laboratory conditions, in accordance with EN 16523-1:2015 (formerly EN 374-3:2003) and EN 374-4:2013. The information is for guidance only and may not reflect the user's application. A risk assessment should always be made by purchaser to assess the suitability of gloves for a specific application.