

## BlotBoy™ Low Speed Gyrotory Rockers (3D)



**Biozym**  
SCIENCE IS OUR BUSINESS



- + Optimized speed & tilt for blotting, washing and staining
- + Tri-directional motion for thorough, gentle mixing
- + Cold room and incubator safe

### Technical Data:

Speed Range:	<b>Fixed, 10 rpm</b>
Tilt Angle:	<b>Fixed, ± 8°</b>
Platform Size (55B3D1008-E):	<b>27 x 19 cm</b>
Platform Size (55B3D1308-E):	<b>30 x 30 cm</b>
Load Capacity:	<b>1.8 kg</b>
Operating Temp. Range:	<b>+4 °C to +65 °C</b>
Dimensions: (W x D x H)	<b>27 x 19 x 20 cm 30 x 30 x 20 cm</b>
Weight:	<b>2 kg</b>
Electrical:	<b>230V, 50-60Hz, 50W</b>
Warranty:	<b>2 Years</b>

Modeled after our popular BioMixer™ series of “nutating” shakers, these platform rockers provide the perfect 3-dimensional motion for all blot related applications. Their tilt angle and speed have been optimally set for gentle, but thorough mixing in gel trays, boxes and other flat vessels.

### Ordering Information:

Item No.	Description:
55B3D1008-E	Mini BlotBoy with flat mat (27 x 19 cm)
55B3D1308-E	BlotBoy with flat mat (30 x 30 cm)
55B3D-DIMPLED	Dimpled mat (27 x 19 cm) for 55B3D1008-E
55BR1000-DIMPLED	Dimpled mat (30 x 30 cm) for 55B3D1308-E
55B0718	Bungee cord, used for securing vessels, pack of 4

## BioMixer™ Nutating Mixer (3D)



- + Includes dimpled mat; perfect for mixing blood & centrifuge tubes
- + Tri-directional motion for thorough, uniform mixing

These multi-functional shakers provide a 3-dimensional “nutating” motion for thorough, foam-free mixing in blood and centrifuge tubes, as well as in a variety of flat vessels.

### Ordering Information:

Item No.	Description:
55B3D1020-E	Mini BioMixer with dimpled mat (27 x 19 cm)
55B3D1320-E	BioMixer with dimpled mat (30 x 30 cm)
55B3D-FLAT	Flat mat (27 x 19 cm) for 55B3D1020-E
55BR1000-FLAT	Flat mat (30 x 30 cm) for 55B3D1320-E
55B0718	Bungee cord, used for securing vessels, pack of 4

### Technical Data:

Speed Range:	<b>Fixed, 20 rpm</b>
Tilt Angle:	<b>Fixed, ± 20°</b>
Platform Size: (B3D1020):	<b>27 x 19 cm</b>
Platform Size: (B3D1320):	<b>30 x 30 cm</b>
Load Capacity:	<b>1.8 kg</b>
Operating Temp. Range:	<b>+4 °C to +65 °C</b>
Dimensions (W x D x H)	<b>27 x 19 x 20 cm 30 x 30 x 20 cm</b>
Weight:	<b>2 kg</b>
Electrical:	<b>230V, 50-60Hz, 50W</b>
Warranty:	<b>2 Years</b>

[www.biozym.com](http://www.biozym.com)