

# Incubation Trays

Convenient, easy to use trays for staining and washing gels and membranes

These specially designed **incubation trays** are ideal for staining and washing electrophoresis gels and membranes.


Choose from three designs, all available in a range of sizes:

- **Opaque** – for light-sensitive applications
- **Transparent** – for easy monitoring of colorimetric staining
- **Traditional** – perfect for blot washes and incubations




## Advantages

- **SMOOTH INTERIOR** – protects membranes and gels from scratches
- **ATTACHED LID** – to protect experiments from dust or debris that can cause speckles on blot images
- **CONVENIENT** – multiple sizes available for different sizes of gels and blots
- **SAVE** – minimize antibody and buffer usage by using appropriately sized trays



Easy visualization of staining



Tray lids protect from dust and debris that cause speckles on Western blots



Save antibody and buffer by using the smallest tray that can hold your blot



Beautiful Western blots every time

## Advansta Solutions

**Staining with a light-sensitive dye or doing a fluorescent Western?**

**Use Advansta's opaque incubation trays**

- Protect your blot or gel from light
- No need to cover your tray with foil
- No need to keep your shaker in a darkened room
- Use for hybridization and washes of fluorescent microscope slides

**Not sure when staining or destaining is complete?**

**Use Advansta's transparent incubation trays**

- Easily observe staining
- Photograph your gel right in the tray

# Incubation Trays

Convenient, easy to use trays for staining and washing gels and membranes



## Ordering Information

Catalog Number	Product		Size
541094clear	Incubation Tray, Transparent body and lid	5x7 cm	5/pack
541094black	Incubation Tray, Opaque black body and lid	5x7 cm	5/pack
541094	Incubation Tray, Traditional red body with clear lid	5x7 cm	5/pack
541095clear	Incubation Tray, Transparent body and lid	6x9 cm	4/pack
541095black	Incubation Tray, Opaque black body and lid	6x9 cm	4/pack
541095	Incubation Tray, Traditional black body with clear lid	6x9 cm	5/pack
541096clear	Incubation Tray, Transparent body and lid	9x11 cm	3/pack
541096black	Incubation Tray, Opaque black body and lid	9x11 cm	3/pack
541096	Incubation Tray, Traditional blue body with clear lid	9x11 cm	5/pack
541097clear	Incubation Tray, Transparent body and lid	10.5x15.5 cm	2/pack
541097black	Incubation Tray, Opaque black body and lid	10.5x15.5 cm	2/pack

## Related Products

541003	WesternBright™ ECL Western Blotting HRP Substrate	Trial kit size	20 ml
541004	WesternBright™ ECL Western Blotting HRP Substrate	For 2000 cm <sup>2</sup> membrane	200 ml
541005	WesternBright™ ECL Western Blotting HRP Substrate	For 5000 cm <sup>2</sup> membrane	500 ml
541089	Goat-anti-mouse HRP-conjugated secondary antibody		500 µl
541088	Goat-anti-rabbit HRP-conjugated secondary antibody		500 µl
541040	Low-Fluorescence PVDF Transfer Membrane	7x9 cm	10 sheets
541041	Nitrocellulose Transfer Membrane 0.45 µm	8x10 cm	10 sheets
541042	Nitrocellulose Transfer Membrane 0.22 µm	8x10 cm	10 sheets

### Advansta Corporation

1505 Adams Drive, Suite B1 | Menlo Park, CA 94025

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: [sales@advansta.com](mailto:sales@advansta.com)

Product information: [www.advansta.com/Incubation\\_trays](http://www.advansta.com/Incubation_trays)

Copyright © 2010–2014 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

[www.advansta.com](http://www.advansta.com)

