

Avant™ buffer pouches

Convenient, pure, fast buffer preparation

Avansta's **Avant buffer pouches** are perfect for standardizing your electrophoresis and Western blotting applications. Dissolve a pouch of pre-measured molecular biology-grade chemicals in 500 ml of deionized water to obtain a ready-to-use 1x buffer without any further preparation. Save precious time in the laboratory and generate reproducible results with Avant buffers.



Advantages

Quality – molecular biology-grade buffers of highest purity

Convenient – pre-measured powder in sealed pouches

Fast – dissolve and use right away

Standardized – guaranteed reproducible results for your electrophoresis gels and Western blots

Ordering Information

Catalog Number	Product	Size
541048	PBS (concentration 1x at 500 ml)	20 pouches
541049	TBS (concentration 1x at 500 ml)	20 pouches
541050	Tris-Glycine PAGE running buffer (concentration 1x at 500 ml)	20 pouches
541051	Tris-Glycine SDS PAGE running buffer (concentration 1x at 500 ml)	20 pouches

Related Products

541030	WesternBright™ MCF fluorescent Western blotting kit (Contains: secondary antibodies, washing solution, blocking solution, pre-cut PVDF membranes and background quenching sheets)	10 blots
541014	WesternBright™ Quantum Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
541004	WesternBright™ ECL Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml

Copyright © 2010 Avansta. All rights reserved. The Avansta logo is a registered trademark of the Company. WesternBright™ and Avant™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Avansta Corporation

1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com

www.advansta.com

