

		Single Tubes				8er und 12er Tube Strips						Strips with attached caps				96-Well Format								Plates 384-Well Format		
Produkt		Art. Nr. 712028	Art. Nr. 712000	Art. Nr. 712010	Art. Nr. 712040	Art. Nr. 712030(8)	Art. Nr. 712098	Art. Nr. 712070	Art. Nr. 712050	Art. Nr. 712150	Art. Nr. 712180	Art. Nr. 712155	Art. Nr. 712120(8)	Art. Nr. 712130	Art. Nr. 712190	Art. Nr. 712195	Art. Nr. 712400	Art. Nr. 712240	Art. Nr. 712200	Art. Nr. 712220	Art. Nr. 712260	Art. Nr. 712420	Art. Nr. 712410	Art. Nr. 712280	Art. Nr. 712450	Art. Nr. 712460
Cyclertyp		RP	RP	HP	LP	RP	RP	HP	LP	LP	LP	RP	RP	HP	LP	RP	RP	LP	RP	RP	LP	RP	RP	LP		
ABI																										
2400	RP	●	●	●		●	●	●				●	●	●		●	●C					●C	●C			
2700 u. 2720	RP	●				●	●	●				●	●	●		●	●			●	●●		●	●		
5700	RP	●				●	●				●	●●	●	●		●	●			●	●		●	●		
7000	RP	●				●	●	●			●	●●	●	●		●	●			●	●		●	●		
7300	RP					●						●●	●			●	●			●●		●●	●●			
7500	RP					●						●●	●			●	●			●●		●●	●●			
7500 Fast	LP				●				●A	●	●A				●A					●●						
7700	RP	●	●			●	●					●●	●			●	●			●●		●	●			
7900	RP											●●				●A				●		●	●		●●	
7900 Fast	LP				●				●A		●A				●A					●●						
9600	RP/HP	●	●	●	●	●	●					●	●		●	●	●			●	●●		●●	●●		
9700	RP	●				●	●					●●	●			●	●			●	●●		●●	●●		●●
9800 Fast	LP				●				●●	●●	●●				●					●●					●	
StepOne	LP				●				●	●	●				●			●C						●C		
StepOne Plus	LP				●				●	●	●				●					●●				●		
Veriti	RP	●				●	●					●●	●			●	●		●	●●		●●	●●		●●	
Veriti Fast	LP				●				●	●	●				●			●		●●						
Biometra																										
Tpersonal	MP 48	●	●	●	●	●	●	●	●	●		●	●	●	●	●	●C	●C				●C	●C			
T1	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●●	●	●			●	●	●	●	●
Tgradient	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●●	●	●			●	●			
T3000	MP	●	●	●	●	●	●	●	●	●		●	●	●	●	●	●C					●C	●C			
Tprofessional	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●●	●	●			●	●		●	●
Trobot	MP	●	●	●	●	●	●	●				●	●	●	●	●	●●	●	●			●	●	●	●	●
BioRad																										
GeneCycler	RP	●	●	●		●	●	●				●●	●	●		●	●C		●			●C	●C			
iCycler	RP	●	●	●		●	●	●				●●	●	●		●	●			●		●	●●		●	●
iCycler IQ	RP	●	●	●●		●	●	●				●●	●	●		●	●			●		●●	●●		●	●
IQ5	RP	●	●	●●		●	●	●				●●	●	●		●	●			●		●●	●●			
MyCycler	RP	●	●	●●		●	●	●				●●	●	●		●	●			●		●●	●●			
MylQ	RP	●	●	●●		●	●	●				●●	●	●		●	●			●		●●	●●			
PTC-100	MP	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●			●	●	●		
PTC-150	MP	●	●	●	●							●	●	●		●	●									
PTC-200 u. 225	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●	●	●
Opticon	LP				●				●●	●●	●●				●			●							●	
Opticon2	LP				●				●●	●●	●●				●			●							●	
MiniOpticon	MP	●	●	●●	●	●	●	●	●●	●●	●●	●●	●	●	●	●	●C					●C	●C			
Chromo4	MP	●	●	●●	●	●	●	●	●●	●●	●●	●●	●	●	●	●	●	●	●			●	●	●		
C-1000	LP				●				●●	●●	●●				●										●	
S-1000	LP				●				●●	●●	●●				●			●							●	
CFX96	LP				●				●●	●●	●●				●			●							●	
CFX384																									●	●
Eppendorf																										
MasterCycler	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●		
MC Gradient	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●		
MC ep Gradient	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●			
MC ep Gradient S	MP	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●		●	●
MC ep Realplex	RP	●	●	●●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●			●	●		●	●
MC Personal 16					●																					
Stratagene																										
Robo-Cycler	MP	●		●	●	●	●	●	●	●	●	●	●	●		●	●	●	●			●	●			
MX3000P	RP	●	●	●●		●	●					●	●			●	●			●		●	●			
MX3005P	RP	●	●	●●		●	●					●	●			●	●			●		●	●			
MX4000		●	●	●		●	●					●				●	●			●		●	●			
Roche																										
LightCycler 480	LP				●A				●A	●A	●A				●A			●A						●●		●●
LightTyper	LP				●A				●A	●A	●A				●A			●A						●●		●●
Analytic Jena / Biozym																										
FlexCycler	RP																●									

red = Low Profile Cycler (96 x 0.1 ml)
LP = Low Profile (0.1 ml Volumen)

RP = Regular Profile
HP = High Profile
MP = Accept all Profiles

• = Passt
•• = Passt ideal
•C = Kann geschnitten werden
•A = Adapter erforderlich

Alle Platten sind auch in „weiß“ erhältlich!
Art. Nr. siehe Tabelle Seite 42.