



96 Well 8-Cap Strip Platte, reißbar 96 Well 0.2 ml 8-Tube Strip Platte, reißbar

FEATURES	ADVANTAGES	BENEFITS
Tear Off capability	Individual closure of tube strips	Minimize cross contamination
	Closure of plates or sections of a plate	Economic usage
	Opening of individual rows of a plate	Minimize cross contamination
Optical clear	Higher S/N ratios qPCR applications	Low Copy detection
Wide area indented caps	No "finger touch" while closing	Highly reproducible results
Extra robust strips	No twisting or bending	Easy to handle
	In line and no breakage	Easy positioning and handling
Clever design cap	Easy opening and closing	"Finger" and RSI Friendly
	Minimal evaporation	Consistent reaction conditions
Auto-orientation holes	Easy sample identification	Traceability samples
Numeric markings	Easy sample location	Traceability samples




User procedures and examples


Tube Strips



Multiple Tube Strips or sections of 96 well plates



96 well plates



Art.Nr.	Bezeichnung	
712380	96 Well 8-Cap Strip Platte, reißbar, farblos, Optical Flat Cap	(Packung à 25 Platten = 300 Deckelstreifen)
High Profile		
712510	96 Well 0.2 ml 8-TubeStrip Platte, reißbar, farblos	(Packung à 25 Stück = 2.400 Tubes)
712511	96 Well 0.2 ml 8-TubeStrip Platte, reißbar, frosted	(Packung à 25 Stück = 2.400 Tubes)
712512	96 Well 0.2 ml 8-TubeStrip Platte, reißbar, weiß	(Packung à 25 Stück = 2.400 Tubes)
Low Profile		
712515	96 Well 0.1 ml 8-TubeStrip Platte, reißbar, farblos	(Packung à 25 Stück = 2.400 Tubes)
712516	96 Well 0.1 ml 8-TubeStrip Platte, reißbar, frosted	(Packung à 25 Stück = 2.400 Tubes)
712517	96 Well 0.1 ml 8-TubeStrip Platte, reißbar, weiß	(Packung à 25 Stück = 2.400 Tubes)