

RNase A

Cat. No. MRNA092



Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook ([facebook.com/EpicentreBio](https://www.facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

Ribonuclease A (RNase A) is an endoribonuclease that cleaves single-stranded RNA at the 3' end of pyrimidine residues, forming oligoribonucleotides having 3'-terminal pyrimidine-3'-phosphates. Pyrimidine-3'-monophosphates are also released by RNase A cleavage of adjacent pyrimidine nucleotides.

Applications

- Removal of RNA from DNA preparations.
- Removal of unhybridized regions of RNA from DNA-RNA or RNA-RNA hybrids.

RNase A is available in a 2 ml size at a concentration of 5 mg/ml.

2. Product Specifications

Storage: Store only at -20°C in a freezer without a defrost cycle.

Composition and Concentration: RNase A is provided on the basis of weight at 5 mg/ml in Storage Buffer.

Storage Buffer: RNase A is supplied in a 50% glycerol solution containing 25 mM sodium acetate (pH 4.6).

Contaminating Activity Assays: RNase A is free of contaminating exo- and endodeoxyribonuclease activities.

3. Related Products

The following products are also available:

- RiboShredder™ RNase Blend
- RNase I, *E. coli*
- RNase III, *E. coli*
- DuraScribe® T7 Transcription Kit

DuraScribe is a registered trademark, and RiboShredder is a trademark of Epicentre, Madison, Wisconsin.

Visit our technical blog: epicentral.blogspot.com