


LucentBlue™ X-ray film

High sensitivity X-ray film for chemiluminescent Western blots

LucentBlue is a sensitive film ideal for imaging chemiluminescent Western blots and is guaranteed to perform with WesternBright™ ECL, WesternBright Quantum and WesternBright Sirius. LucentBlue is compatible with horseradish peroxidase and alkaline phosphatase detection reagents, as well as autoradiography. In addition to Western blots, LucentBlue can be used for Northern and Southern blot detection, and other techniques requiring film.



Advantages

- **Performance** – excellent images when used with WesternBright ECL, WesternBright Quantum, and WesternBright Sirius
- **Value** – sensitivity equal to more expensive films
- **Flexibility** – detect radioisotopes in addition to chemiluminescence

Ordering Information

Catalog Number	Product	Size
() (%\$- \$.....)	LucentBlue™ X-ray film, 8x10" sheets	100 sheets
() (%\$- %.....)	LucentBlue™ X-ray film, 5x7" sheets	100 sheets

Related Products

541004	WesternBright™ ECL Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541015	WesternBright™ Quantum Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541021	WesternBright™ Sirius Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541089	Goat-anti-mouse HRP-conjugated secondary antibody	500 µl
541088	Goat-anti-rabbit HRP-conjugated secondary antibody	500 µl
541040	Low-Fluorescence PVDF Transfer Membrane 7x9 cm	10 sheets
541041	Nitrocellulose Transfer Membrane 0.45 µm 8x10 cm	10 sheets
541042	Nitrocellulose Transfer Membrane 0.22 µm 8x10 cm	10 sheets

Copyright © 2010 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ and LucentBlue™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Advansta Corporation
1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025
Tel.: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com
www.advansta.com

