

Accuris™ Mini Balances



- + Super compact, less than 7 cm wide
- + Easy to read, backlit LCD
- + Calibration function
- + Includes two AAA batteries
- + Snap on, protective cover



Optional in-lab
calibration weights



Protective cover,
included

Accuris Mini Balances are extraordinarily compact, easy to use, and convenient for a variety of laboratory weighing applications. With a footprint of only 6.6 x 9.9 cm, these Mini Balances offer exceptional accuracy and linearity through their full weighing ranges. Available with either a 100g or 500g capacity, these digital balances feature an easy to read, backlit, liquid crystal display. Buttons for

ON/OFF, MODE & TARE are located on the front for easy access. The software includes a calibration function for in-lab calibration with a 100 g or 500 g calibration weight (available separately).

Accuris Mini Balances are powered by 2 AAA batteries (included). After one minute sitting idle, the power shuts off to conserve battery life. Ideal for lab or field use, each balance includes a snap on, protective cover with instructions conveniently printed on the underside.

Model	Mini 100 g	Mini 500 g
Max capacity	100 grams	500 grams
Readability	0.01 g	0.1 g
Calibration	External	
Pan size	5.3 x 5.7 cm	
Weighing units	g, ct, dwt, gn	
Display	backlit LCD	
Power	2 AAA batteries	
Exterior dims.	6.6 x 9.9 x 1.9 cm	

Ordering Information:

Item No.	Description:
55W4000-100	Accuris™ Mini Balance, 100 g capacity (0.01 g readability)
55W4000-500	Accuris Mini Balance, 500 g capacity (0.1 g readability)

Accuris™ Calibration Weights



- + Various individual and weight sets available
- + Choose E2, F1, or M1 OIML class depending on requirements
- + Weight sets supplied in convenient protective cases

Accuris calibration weights are ideal for routine calibration of your Accuris (or any brand) lab balance. Single weights as well as complete weight sets are available in various tolerance classes: OIML class E2, F1, and M1.

M1 Class: Ideal for basic (compact) 0.1 to 0.01 g readability balances
F1 Class: Suitable for precision 0.01 to 0.001g readability balances
E2 Class: Recommended for analytical 0.0001g readability balances

Individual weights come in a protective case. Weight sets include a protective carrying case with foam liner.

www.biozym.com

Ordering Information:

55W1000-100	Accuris™ Calibration Weight (F1 grade) 100 g
55W1000-200	Accuris™ Calibration Weight (F1 grade) 200 g
55W1000-500	Accuris™ Calibration Weight (F1 grade) 500 g
55W1005-1000	Accuris™ Calibration Weight (M1 grade) 1 kg
55W1005-2000	Accuris™ Calibration Weight (M1 grade) 2 kg
55W1000-TF	Tweezers, for handling calibration weights
55W1101-6-100	Calibration set for Analytical Balances (Class E2) 6 pcs: 1x100 mg, 1x10 g, 2x20 g, 1x50 g, 1x100 g, incl. protective case, tweezers & cleaning cloth
55W1100-200	Calibration set for Precision Balances (Class F1) 23 pcs, 1mg to 200g, includes protective case, tweezers, dust brush and cleaning cloth
55W1105-9-1000	Economy calibration set for compact balances (Class M1), 9 pcs: 1 x 10 g, 2 x 20 g, 1 x 50 g, 1 x 100 g, 2 x 200 g, 1 x 500 g, 1 x 1 kg, includes protective case