

TubeRocker™ Fixed Speed Rocker (2D)



Biozym
SCIENCE IS OUR BUSINESS



CE

- Ideal for tube sizes 1.5 ml to 50 ml
- Grooved mat properly secures tubes without the need for any tools, clamps or hardware

Technical Data:


Speed Range:	Fixed, 24 rpm (230V: 20 rpm)
Tilt Angle:	±30°
Platform Size:	30 x 11.4 cm
Load Capacity:	2.3 kg
Operating Temp. Range:	+4 °C to +45 °C
Dimensions: (W x D x H)	33.6 x 13.5 x 10 cm
Weight:	1.8 kg
Electrical:	230V, 50-60Hz, 50W
Warranty:	2 Years



8x50ml



15x15ml



15 x Blood Collection Tubes
(5ml, 7ml, 10ml or 15ml)



30 x 1.5/2.0ml

Ordering Information:

Item No.	Description:
55M2100-E	TubeRocker with grooved mat

BenchBlotter™ Fixed Speed Rocker (2D)



CE



55BR1000-E with
55BR1000-STACK

- Gentle mixing for blot washing and gel staining
- Optional dual platform doubles work surface

Technical Data:

Speed Range:	Fixed, 12 rpm Fixed, 10 rpm (230V)
Tilt Angle:	±8°
Platform Size:	30 x 30 cm
Load Capacity:	2.3 kg
Operating Temp. Range:	+4 °C to +45 °C
Dimensions (W x D x H):	34 x 30 x 12 cm
Stacking Clearance:	7.6 cm
Weight:	3 kg
Electrical:	230V, 50-60Hz, 50W
Warranty:	2 Years

Ordering Information:

Item No.	Description:
55BR1000-E	BenchBlotter Platform Rocker with flat mat
55BR1000-STACK	Optional stacking platform with flat mat
55BR1000-DIMPLED	Optional dimpled mat
55B0718	Bungee cord, used for securing vessels, pack of 4