

Protein Sample Loading Buffers

Pre-mixed loading buffers for polyacrylamide electrophoresis

Advansta provides ready-to-use **Laemmli loading buffers** for the preparation of protein samples for SDS-polyacrylamide gel electrophoresis. Choose from reducing or non-reducing buffers; each contains Bromophenol blue, to enable visualization of the progress of electrophoresis by observing migration of the dye front. Both loading buffers are compatible with downstream Coomassie and silver staining and Western blotting applications.



Advantages

- **Convenient** – Pre-mixed; just add an equal volume to your protein sample, boil and load
- **Quality controlled** – Obtain reproducible results with standardized buffers
- **Flexibility** – Choose from reducing or non-reducing loading buffers

Ordering Information

Catalog Number	Product	Size
541100	Non-reducing protein sample loading buffer (2x)	1 ml
541101	Non-reducing protein sample loading buffer (2x)	5 ml
541102	Reducing protein sample loading buffer (2x)	1 ml
541103	Reducing protein sample loading buffer (2x)	5 ml

Related Products

541030	WesternBright™ MCF fluorescent Western blotting kit (Contains: secondary antibodies, washing solution, blocking solution, pre-cut PVDF membranes and background quenching sheets)	10 blots
541014	WesternBright™ Quantum Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
541020	WesternBright™ Sirius Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
541004	WesternBright™ ECL Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml

Copyright © 2011 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Advansta Corporation
1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025
Tel.: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com
www.advansta.com

