

ORBI-BLOTTER™



- + Gentle orbital motion for gel applications
- + Large 35 x 30 cm work surface
- + Incubator and cold room safe
- + Variable speed control, 3 to 70 rpm

The Orbi-Blotter™ features a gentle orbital motion with speed control down to as little as 3 rpm, making it the ideal orbital choice for low speed applications, such as blotting, washing, staining/destaining, etc. A large 35 x 30 cm work surface is included along with a non-slip rubber mat for holding trays, dishes, plates and other flat vessels.

An optional stacking platform can be purchased as an accessory, effectively doubling the usable work surface. The maintenance free, brushless motor is capable of loads up to 4 kg and is safe for use in cold rooms and incubators.

Ordering Information:

Item No.	Description:
55BT30-E	Orbi-Blotter™ Low Speed Orbital Shaker
55BR2000-STACK	Stacking platform with 35 x 30 cm flat mat (7.6 cm clearance)
55BR2000-SP	STACK-PLUS Extender set, adds 3.2 cm clearance to stacking platform

Technical Data:

Speed Range:	Variable, 3 to 70 rpm
Motion / Orbit:	Orbital / 19 mm
Load Capacity:	3 x 2 l or 4 kg
Platform Size:	35 x 30 cm
Operating Temp. Range:	+4 °C to +45 °C
Dimensions: (W x D x H)	35.5 x 30.5 x 20 cm
Stacking Clearance:	7.6 cm
Weight:	3.4 kg
Electrical:	230V, 50-60Hz, 50W
Warranty:	2 Years

See Pg. 44-45

Shaking Incubators



INCUSHAKER™ Series

- + Exceptional temperature uniformity
- + Three models,
 - Incu-Shaker Mini: Up to 5 x 500 ml (or 2 x 1 l)
 - Incu-Shaker 10 L: Up to 5 x 2 l
 - Incu-Shaker 10 LR: Up to 5 x 2 l (with cooling)
- + Exchange flask clamps in one second with **MAGic Clamp™** magnetic system