



PRACTICAL INNOVATION

Biozym
SCIENCE IS OUR BUSINESS

Ziath

DataPaq™ Linear Rack Scanner Version 2

2-year guarantee as standard

The Ziath Linear Rack Scanner Version 2 is for use with the DataPaq™ Express, Deep Focus, High Speed and Multi-Rack Scanners. It is simply attached to the scanner using a bracket with a magnetic place holder. No assembly of the linear scanner is required.



Features

- Use with Ziath DataPaq™ 2D barcode rack scanners
- Very quick and easy to install. No assembly required.
- Scans in milliseconds
- Wide range of placement options
- Quick & easy installation, no dongles required
- USB connection
- Operating systems: Windows XP, Vista, 7
- Attaches to the rack scanner using bracket with magnetic place holder
- Dimensions : 67 x 48 x 19mm (WxHxD)
- Durable construction - 2 year guarantee
- Free software upgrades for life
- Ziath's expertise and support come as standard



DataPaq™ Linear Rack Scanner Version 2

PRACTICAL INNOVATION

Easy Set-Up and Easy To Use

The DataPaq™ Linear Rack Barcode Scanner is very simple to attach to existing 2D scanners, and offers easy installation. It can be attached to any side of the DataPaq™ scanners using a bracket with a magnetic pad.

Our DataPaq™ software makes it simple to export the linear (rack) barcode data and the

tube data to Excel, XML, JSON or text and scanned images can also be saved.

The linear rack scan takes less than one second and the scanner is powered and controlled by USB cable, so needs no additional power cables.

The Deep Focus Linear Rack Scanner installed on the Deep Focus 2D barcode scanner



Linear Rack Scanner shown with the DataPaq™ Express Rack Reader (top) and the High Speed Rack Scanner (bottom)

Ordering Information

Code	Description	Quantity
601122	ZTS-1DR2-A6 DataPaq™ High Speed Linear Rack Scanner Version 2	1
601121	ZTS-1DR2-AV6 DataPaq™ Deep Focus Linear Rack Scanner Version 2	1
601123	ZTS-1DR2-EXP DataPaq™ Express Linear Rack Scanner Version 2	1
	ZTS-1DR2-A4 DataPaq™ Multi-Rack Linear Rack Scanner Version 2	1

2012_1



Ziath Ltd, The Studio, High Green, Great Shelford, Cambridge, CB22 5EG, UK

Tel: +44 (0) 1223 655577 Fax: +44 (0) 1223 790069 Email: info@ziath.com

www.ziath.com



Biozym Scientific GmbH Tel.: +49 (0)5152-9020 - Fax: +49 (0)5152-2070 support@biozym.com www.biozym.com