



Product Information Sheet N2050-10-L Rev A

Product Information	
10mM dNTP Solution Mix	
Part Number	N2050-10-L 280600
Concentration	10 mM
Unit Size	100 µmol

Product Specification	
Storage Temperature	-25°C to -15°C
TEST:	SPECIFICATION:
Purity	≥99%
Base Purity	≥99.5%
Pyrophosphate	≤0.003 pmol/pmol of nucleotide
ATP	None detected
Nick Activity	None detected
DNase Activity	None detected
RNase Activity	None detected

Product Description:

An equimolar (10 mM) mixture of the four natural deoxynucleotides.

Supplied as:

10 mM each dNTP (dATP, dCTP, dGTP, dTTP) at pH 7.5



Biozym Scientific GmbH
 Steinbrinksweg 27
 D-31840 Hessisch Oldendorf
 Tel.: +49 (0)5152-9020
 Fax: +49 (0)5152-2070
 support@biozym.com
 www.biozym.com



Limitations of Use

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.