


**Green Elephant Biotech GmbH**  
 Kerkrader Str. 9  
 35394 Gießen  
 Germany  
[info@greenebt.com](mailto:info@greenebt.com)


**Tolerances:**  
 Either depicted in the  
 corresponding sketch  
 itself or in the table  
 below.

**Green Elephant®**  
**96-well microPLAtes**

Sheet 1/1  
 31.10.2025



Tolerances:  
 $\pm 0.1$  mm



1. Dimensions are in Millimeters
2. Refer to additional product documents for further information.

Dimension without tolerances (mm)		
0-10	$\pm 0.1$	based on DIN 16901 140B
10-22	$\pm 0.17$	
22-40	$\pm 0.24$	
40-70	$\pm 0.33$	
70-90	$\pm 0.5$	