


AdvanWash™ Washing Solution

Washing solution for chemiluminescent and fluorescent Western blots

AdvanWash washing solution can be used with any Western blotting system, including all WesternBright chemiluminescent HRP substrates and the WesternBright MCF fluorescent Western blotting kit. Provided as a 10X concentrate.

Advantages

- **Reliable** – Quality-controlled wash buffer for reproducible results
- **Versatile** – Compatible with chemiluminescent and fluorescent Western blots
- **Convenient** – Save time using pre-made buffers. No more weighing, mixing, or pH adjusting.



Ordering Information

Catalog Number	Product	Size
541004	AdvanWash™ Washing Solution	500 ml

Related Products

541030	WesternBright™ MCF fluorescent Western blotting kit (Contains: secondary antibodies, washing solution, blocking solution, pre-cut PVDF membranes and background quenching sheets)	10 blots
541014	WesternBright™ Quantum Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
541015	WesternBright™ Quantum Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541020	WesternBright™ Sirius Western Blotting HRP Substrate (for 1000 cm ² membrane)	100 ml
541021	WesternBright™ Sirius Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541004	WesternBright™ ECL Western Blotting HRP Substrate (for 2000 cm ² membrane)	200 ml
541005	WesternBright™ ECL Western Blotting HRP Substrate (for 5000 cm ² membrane)	500 ml
541006	WesternBright™ ECL Spray	500 ml

Copyright © 2011 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. AdvanWash™ and WesternBright™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Advansta Corporation
 1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025
 Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com
www.advansta.com

