



Ziath Handheld 2D Barcode Scanner



Features

- Intelligent and fully portable 2D barcode reader for single tubes
- 4 operating modes
- Up or Down-load data from host computer
- Full day operation without the need to re-charge the battery
- Cryoprotection provided as standard

Four Operating Modes



- 1) Simple barcode scanning**
Display the barcode value from a scanned tube
- 2) Check sample**
Display complete sample details from a pre-loaded check list
- 3) Pick List**
Download a list of samples to the Handheld and check them off one at a time
- 4) Create Sample**
Add a new sample based on descriptors downloaded from the host computer

Included in the Handheld Package

- Handheld instrument
- Battery charger
- USB lead to connect instrument to host computer
- Software package to provide user friendly interface between computer and instrument

Specifications

- Dimensions = 155mm x 65mm x 38mm (length x width x max depth)
- Weight = 150g
- Battery charger input = 100-240VAC 60/50Hz, 260mA
- Battery charger output = 5VDC, 2.4A
- Battery life = Typically 8 hours between charges
- Data storage on instrument = Typically up to 60 000 samples
- Display size = 55mm x 90mm
- Instrument operating system = Android 5.0 Lollipop

Ordering Information

Code	Description	Quantity
601090	Ziath Intelligent Handheld 2D barcode reader	1



Biozym
SCIENCE IS OUR BUSINESS

Biozym Scientific GmbH
Steinbrinksweg 27
D-31840 Hessisch Oldendorf
Tel.: +49 (0)5152-9020
Fax: +49 (0)5152-2070
support@biozym.com
www.biozym.com