



CoolCube™

- + Maintains 0 °C for 4 hours
- + Ideal for:
 - 1.5 ml & 2.0 ml microtubes
 - 0.2 ml PCR tubes & strips
 - 96 well PCR plates
- + Also for use as a PCR workstation

Technical Data:

Temperature:	0 °C for up to 4 hours
Capacity:	
Top Side:	36 x 1.5/2.0 ml tubes
Bottom Side:	96 x 0.2 ml tubes 12 x 0.2 ml strips 1 x 96 well PCR plate 10 x 1.5/2.0 ml
Dimensions: (W x D x H)	14.3 x 14.3 x 6.4 cm
Weight:	1 kg
Warranty:	2 Years

The CoolCube™ is designed to keep samples cold on the lab bench without the mess of a melting ice bucket. Simply place the CoolCube in a freezer overnight and it's ready for use, keeping samples cold, but not frozen. Maintains temperature at approximately 0 °C for up to four hours.

Ordering Information:

Item No.	Description:
55R1000	CoolCube Microtube and PCR Plate Cooler

hybex™ Media Storage Bottles



hybex™

- + Borosilicate 3.3 hybex™ glass
- + Autoclavable caps (GL45) & drip rings
- + Permanent enamel graduation markings
- + Large area for bottle marking/identification
- + Autoclavable, reuseable & breakage resistant
- + Five cap colors

Technical Data:

Material:	hybex borosilicate glass (3.3)
Cap Thread Size:	GL45
Temperature Range:	-70° to 500 °C
Sterilization:	Autoclave
Cap Material:	Polypropylene (140 °C)
Cap Colors:	Standard blue, green, purple, yellow, white or assorted

Ordering Information:

Item No.	Description:
55B3000-50-B	Bottle, 50 ml with GL32 blue cap, 10/pk.
55B3000-100-B*	Bottle, 100 ml with standard (GL45) blue cap, 10/pk.
55B3000-250-B*	Bottle, 250 ml with standard (GL45) blue cap, 10/pk.
55B3000-500-B*	Bottle, 500 ml with standard (GL45) blue cap, 10/pk.
55B3000-1000-B*	Bottle, 1 l with standard (GL45) blue cap, 10/pk.
55B3000-2000-B*	Bottle, 2 l with standard (GL45) blue cap, 10/pk.
55B3000-SP-B*	Bottle, Starter Pack with standard (GL45) blue cap, includes, 2 x 100 ml, 3 x 250 ml, 3 x 500 ml and 2 x 1 l)
55B3000-CAP-B*	Replacement cap, standard blue (GL45), 10/pk.
55B3000-CAP2	Replacement cap, standard blue (GL32), 10/pk.
55B3000-CAP-HTC	High temperature caps, red (GL45), up to 180°, 10/pk.