



55H1000-MR-CMB Adjustable flask clamp, accepts 125 ml, 250 ml and 500 ml Erlenmeyer flasks.

(For a video demonstration, scan the QR code above.)

Maximum Flask Capacity with MAGic Clamp Platforms

Item No.	Description:	Platform 55H1000-MR	Platform 55BT3000-MR	Platform 55H1010-MR
		Available For Orbi-Shaker JR.	Available For Models Orbi-Shaker & Orbi-Shaker CO ₂	Available For Orbi-Shaker XL (CO ₂)
55H1000-MR-CMB	NEW! Adjustable flask clamp, accepts 125, 250 & 500 ml Erlenmeyer flasks	12 / 8 / 5	20 / 14 / 9	36 / 25 / 16
55H1000-MR-50	Flask clamp, 50 ml Erlenmeyer flask	20	30	49
55H1000-MR-125	Flask clamp, 125 ml Erlenmeyer flask	12	20	36
55H1000-MR-250	Flask clamp, 250 ml Erlenmeyer flask	8	14	25
55H1000-MR-500	Flask clamp, 500 ml Erlenmeyer flask	5	9	16
55H1000-MR-1000	Flask clamp, 1 l Erlenmeyer flask	4	5	9
55H1000-MR-2000	Flask clamp, 2 l Erlenmeyer flask	1	2	5
55H1000-MR-3000	Flask clamp, 3 l (or 2.8 l) Fernbach flask	N/A	1	4

Maximum Capacity for MAGic Clamp Platforms (Test Tube / Microplate Racks)

Item No.	Description:	Platform 55H1000-MR	Platform 55BT3000-MR	Platform 55H1010-MR
		Available For Orbi-Shaker JR.	Available For Models Orbi-Shaker & Orbi-Shaker CO ₂	Available For Orbi-Shaker XL (CO ₂)
55H1000-MR-T15	Tube Rack, 15 ml conical tubes, 32 x 15 ml	2	3	8
55H1000-MR-T50	Tube Rack, 50 ml conical tubes, 15 x 50 ml	2	3	8
55H1000-MR-1550	Combi-Tube Rack, 30 x 15 ml & 20 x 50 ml tubes	2	2	4
55H1000-MR-T600	Tube/Bottle Rack, 4 x 500/600 ml conical bottles	1	2	4
55H1000-MR-T15H	Tube rack, 12 x 1.5/2.0 ml tubes, horizontal	2	4	6
55H1000-MR-T150H	Tube rack, 4 x 15 ml tubes, horizontal	4	6	8
55H1000-MR-T500H	Tube rack, 2 x 50 ml tubes, horizontal	4	6	8
55H1000-MR-MP	Plate rack, holds one microplate	4	6	8
55H1000-MR-TSR	Angled tube rack holder (with adjustable tilt)	2	6	8