

# High-Capacity Sample Transport Tote

## DuraPorter® XL High-Capacity Sample Transport Tote

Your reliable solution for the safe and convenient transportation of laboratory samples, instruments, and valuable items. Crafted to meet the stringent demands of clinical and research laboratories, this oversized transport container ensures the protection and preservation of your cargo.

### Built to Last and Designed with the User in Mind

- Sturdy and spacious, ideal for transporting a large volume of samples
- Transparent design helps comply with OSHA 29 CFR 1910.1030 standards without having to re-label containers and enabling the user to visibly check sample integrity and meets UN3373 secondary packaging requirements
- Equipped with articulated handles featuring silicone grips for comfortable and secure handling



### Emphasizing Security and Safety

- A silicone seal ensures a tight enclosure, safeguarding samples from most environmental elements
- Dual latches with security tag feature guarantee a secure closure, providing tamper-evidence and ensuring the integrity and authenticity of your samples



### Exceptional Versatility

- Includes a MultiTube Rack Insert to accommodate various tube sizes. *Additionally, two other insert options for blood tubes up to 10 mL and general research tubes up to 15 mL are available separately*
- The bottom compartment is generously sized, allowing for carrying most tube racks, instruments or supplies according to your specific requirements

**Soft Grip Handle**  
for comfortable handling

**Translucent Design**  
allows for easy content  
identification

**Removable Lid**  
for ease of loading and  
unloading

**4-Point Handle Security**  
ensures durable transport

**Silicone Seal**  
for secure content enclosure




*Image shown with Small Tube Rack Insert (not included)*

**Dual Latches**  
with tag out feature for  
added security

**Available in 3 Colors**  
for easy user or sample  
identification and organization



Includes One 762655 MultiTube Rack Insert

Item No.	Color		Dimensions (not including handle) cm
761505	Blue		40.64 x 41.28 x 17.78
761506	Red		40.64 x 41.28 x 17.78
761507	Yellow		40.64 x 41.28 x 17.78

DuraPorter® XL  
MultiTube Rack Insert



- Designed to fit DuraPorter® XL High-Capacity Sample Transport Tote
- Thermoformed PETG
- Tube Rack Inserts can be used separately
- Lid removes completely for easy cleaning
- Insert deep enough to be filled with ice

Item No.	Tubes	Cut-out Diameter	UOM
<b>762656</b>	15 mL, Large blood tube	18 mm	Pack of 3
<b>762657</b>	Small blood tubes	13 mm	Pack of 3
<b>762655</b>	50 mL, 15 mL, Small & Large blood tubes	30 mm (center), 18 mm (large), 13 mm (small)	Pack of 3



**762656**

DuraPorter® XL Large Tube Insert

Holds 100 tubes with diameter up to 18 mm  
Most 10 mL blood tubes & 15 mL conical tubes  
Array 10 x 10  
Pack of 3



**762657**

DuraPorter® XL Small Tube Insert

Holds 168 tubes with diameter up to 13 mm  
Most 3.5 - 7.5 mL blood tubes  
Array 12 x 14  
Pack of 3



**762655**

DuraPorter® XL MultiTube Rack Insert

Holds tube with various diameters: 13, 18 & 30 mm  
Most 15 mL & 50 mL conical tubes  
and 3.5 mL - 10 mL blood tubes  
Pack of 3