

---

<b>Product Name:</b>	<b>Green Elephant® 96-well microPLAte</b>	clear	white	black
<b>Product Number:</b>	GEB-WP96FBT	GEB-WP96FBW	GEB-WP96FBB	

---

Green Elephant Biotech certifies that the product has been manufactured in accordance with established manufacturing guidelines and product specifications. The product meets all Green Elephant Biotech quality requirements.

**ISO Certification** - Green Elephant Biotech's quality management system is certified according to ISO 9001:2015-09.

**Functional Testing** - Products undergo continuous quality controls regarding function, precision, and reproducibility. Dimensional checks, testing of precision and accuracy, transparency, etc. are part of Green Elephant Biotech's quality assurance standards.

**Biocompatibility Testing** - Biological testing procedures ensure that the product is free from endotoxins (<0.0125 EU/ml) and does not have a cytotoxic effect in vitro.

Testing is conducted by an independent, EN ISO/ IEC 17025-accredited laboratory. Endotoxin is tested according to EN ISO 10993-11 (2018) and Eu.-Phr.2.6.14. In vitro cytotoxicity is tested according to EN ISO 10993-5 (2009).

**DNA, DNase and RNase Testing** - The product has been tested and confirmed to be free from

- human DNA (Realtime-PCR, LOD: < 2 pg/µl)
- DNase (DNA Digestion, LOD: < 10<sup>-6</sup> Kunitz-Units) and
- RNase (RNA Digestion, LOD: < 10<sup>-9</sup> Kunitz-Units).

Testing is conducted by an independent, EN ISO/ IEC 17025-accredited laboratory.

**Traceability** - The product is fully traceable by lot number.

Sabrina Seidler  
Head of Quality Management