



- + The only vortexer available with **Multi-Head™**
- + Powerful motor for **INSTANT** vortexing
- + **Q-Drive™** dynamic balancing system
- + Variable speed from 200 to 3,200 rpm




Unique Multi-Head allows for standard vortexing, **PLUS** holds up to 8 microtubes

The Mortexer™ is MORE than your standard vortex mixer. The unique design of the Multi-Head™ (included) features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes. The drive system is balanced with Q-Drive™ technology, providing a smooth vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation.

Other features include: "no-walk" suction feet, cold room/incubator compatibility and spill-proof electronic controls.

Available Accessories:



Technical Data:

Speed Range:	200 to 3,200 rpm
Operating Modes:	Touch or continuous
Orbit:	3 mm
Operating Temp. Range:	+4 °C to +45 °C
Dimensions: (W x D x H)	13 x 16 x 17 cm
Weight:	3.8 kg
Electrical:	230V, 50Hz*, 180W
Warranty:	2 Years

*230V, 60Hz available upon request

Ordering Information:

Item No.	Description:
55BV1005-E	Mortexer Vortex Mixer (with Multi-Head)
55BV1000-FLAT	Optional 7.6 cm flat head attachment
55BV1000-COMBO	Optional combination head for tubes and microplates (38x1.5/2.0 ml & 28 x 0.5 ml or one microplate)
55BV1000-H15	Optional head for 12 x 1.5/2.0 ml, horizontal
55BV1000-H150	Optional head for 4 x 15 ml tubes, horizontal
55BV1000-H500	Optional head for 2 x 50 ml tubes, horizontal