

HighPrep DTR



50 ml
Artikel-Nr.: 220022 | MagBio

1.008,90 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Jubiläumsaktion
Gewebe/Organismus: Sequenzierungsansatz
Methode: Magnetische Beads
Nukleinsäure: DNA
Verpackung: 50 ml

Vorteile:

- Effektives Entfernen von Kontaminationen in Sequenzierungsansatz
- Long Phred20 Reads Leselänge >700 Basen
- Hohe durchschnittliche Signalstärke
- Reduzierte Menge BigDye™ möglich

*BigDye is a registered trademark of Applied Biosystems

Anwendung:

Aufreinigung von Sequenzierungsansätzen für ABI und MegaBACE Plattformen

Vergleichende Analyse

HighPrep™ DTR vs CleanSEQ®

Platform	HighPrep DTR	CleanSEQ
Phased DNA	100%	100%
PCR: 1,5kb PCR	100%	100%
Phased DNA	100%	100%
PCR: 1,5kb PCR	100%	100%
Phased DNA	100%	100%
PCR: 1,5kb PCR	100%	100%

Packaging

Volume Size	50 mL
-------------	-------

Packaging	
Storage and Stability	
Storage	Store at 4°C DO NOT FREEZE
Stability	Stable for 12 months when
Disclaimers	
Use Restriction	For Research Use Only. Not for Use in Diagnostic Procedures.
Number of Reactions - 96 Well Format	
5 µL Reaction Volume	5,000
10 µL Reaction Volume	5,000
15 µL Reaction Volume	5,000
20 µL Reaction Volume	5,000
Number of Reactions - 384 Well Format	
5 µL Reaction Volume	10,000
10 µL Reaction Volume	10,000
15 µL Reaction Volume	10,000

CleanSEQ is a registered trademark of Beckman Coulter