

Purified Mouse IgG2b, kappa Isotype Ctrl

Mouse IgG2b, kappa
MPC-11

BioLegend
From Revvity



50 µg

Artikel-Nr.: B400301 | **BioLegend** | Herstellernummer: 400301

24,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ChIP, FA, FC, ICC, ICFC, IHC, IP, WB

Clone: MPC-11

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 50 µg

Purified Mouse IgG2b, kappa Isotype Ctrl

Clone: MPC-11

Isotype: Mouse IgG2b, kappa

Other_Names:

Size: 50 µg

Reactivity:

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA, ChIP

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1416>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>