

## Purified anti-human CD47

Mouse IgG1, kappa  
CC2C6

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B323102 | **BioLegend** | Herstellernummer: 323102

**95,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FC, ICC

**Clone:** CC2C6

**Format:** PURE

**Isotype:** kappa, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 100 µg

Purified anti-human CD47

Clone: CC2C6

Isotype: Mouse IgG1, kappa

Other\_Names: Rh-associated protein, gp42, integrin-associated protein, IAP, neurophilin

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=3705>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>