

PE anti-human MICA/MICB

Mouse IgG2a, kappa
6D4

BioLegend
From Revvity



100 tests

Artikel-Nr.: B320906 | **BioLegend** | Herstellernummer: 320906

190,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: 6D4

Format: PE

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 100 Tests

PE anti-human MICA/MICB

Clone: 6D4

Isotype: Mouse IgG2a, kappa

Other_Names: PERB11, MIC, MHC class I chain-related protein A, MICA, PERB11.1, MHC class I chain-related protein B, MICB, PERB11.2

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=3064>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>