

Purified anti-human CD266 (Fn14, TWEAK Recepto...

Mouse IgG1, kappa
ITEM-1

BioLegend
From Revvity



100 µg

Artikel-Nr.: B314002 | **BioLegend** | Herstellernummer: 314002

138,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: Apop, FC, IHC-P, WB
Clone: ITEM-1
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD266 (Fn14, TWEAK Receptor)

Clone: ITEM-1
Isotype: Mouse IgG1, kappa
Other_Names: Fibroblast growth factor-inducible 14 kD protein (Fn14), TNFRSF12A
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, Apop, WB, IHC-P
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=2665>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>