

Purified anti-human CD255 (TWEAK)

Mouse IgG3, kappa
CARL-1

BioLegend
From Revvity



100 µg

Artikel-Nr.: B308302 | **BioLegend** | Herstellernummer: 308302

117,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC, IHC

Clone: CARL-1

Format: PURE

Isotype: kappa, Mouse IgG3

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human CD255 (TWEAK)

Clone: CARL-1

Isotype: Mouse IgG3, kappa

Other_Names: TNF-like weak inducer of apoptosis, Apo-3 ligand (Apo3L), Death receptor 3-L (DR3-L), TNFSF 12, CD255

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, IHC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1333>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>