

Purified anti-mouse TER-119/Erythroid Cells

Rat IgG2b, kappa

TER-119

BioLegend
From Revvity



50 µg

Artikel-Nr.: B116201 | **BioLegend** | Herstellernummer: 116201

34,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, IHC-F, IHC-P, IP, WB

Clone: TER-119

Format: PURE

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 50 µg

Purified anti-mouse TER-119/Erythroid Cells

Clone: TER-119

Isotype: Rat IgG2b, kappa

Other_Names: Ly-76

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, IHC-F, IHC-P, IP, WB

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1869>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>