

## Purified anti-TCF1 (TCF7)

Rat IgG2a, kappa  
W16175A

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B615701 | **BioLegend** | Herstellernummer: 615701

**81,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** IP, WB

**Clone:** W16175A

**Format:** PURE

**Isotype:** kappa, Rat IgG2a

**Species of Origin:** Rat

**Verpackung:** 25 µg

Purified anti-TCF1 (TCF7)

Clone: W16175A

Isotype: Rat IgG2a, kappa

Other\_Names: Transcription Factor 7, Transcription Factor 7 (T-Cell Specific, HMG-Box), T-Cell-Specific Transcription Factor 1, T-Cell Factor-7, T-Cell Factor 1, TCF1

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App\_Abbrev: WB, IP

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=18183>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>