

Purified anti-Acetylated Lysine

Mouse IgG2b, kappa
15G10

BioLegend
From Revvity



25 µg

Artikel-Nr.: B623401 | **BioLegend** | Herstellernummer: 623401

60,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICC, WB

Clone: 15G10

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-Acetylated Lysine

Clone: 15G10

Isotype: Mouse IgG2b, kappa

Other_Names: Proteins are acetylated by the reversible transfer of acetyl-CoA to the epsilon-amino group of lysine

Size: 25 µg

Reactivity: Species independent

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=2509>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>