

## Purified anti-Eg5

Mouse IgG1, kappa  
10C7/Eg5

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B627802 | BioLegend | Herstellernummer: 627802

**184,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** ICC, IP, WB

**Clone:** 10C7/Eg5

**Format:** PURE

**Isotype:** kappa, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 100 µg

Purified anti-Eg5

Clone: 10C7/Eg5

Isotype: Mouse IgG1, kappa

Other\_Names: Kinesin-like protein KIF11, kinesin-related motor protein Eg5, kinesin-like spindle protein HKSP, thyroid receptor interacting protein 5 (TRIP5)

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App\_Abbrev: WB, ICC, IP

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=2930>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>