

## Purified anti-AND-1 (WDHD1)

Mouse IgG2a, kappa  
20G10

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B630301 | [BioLegend](#) | Herstellernummer: 630301

**95,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** IHC-P, WB

**Clone:** 20G10

**Format:** PURE

**Isotype:** kappa, Mouse IgG2a

**Species of Origin:** Mouse

**Verpackung:** 25 µg

Purified anti-AND-1 (WDHD1)

Clone: 20G10

Isotype: Mouse IgG2a, kappa

Other\_Names: WD repeat and HMG box DNA binding protein 1, acidic nucleoplasmic DNA-binding protein 1

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App\_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=3922>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>