

PCR TW 8-Tube Strip, 0.2 ml, frosted

High Profile, Extra Robust

DNase-, RNasefrei



Beutel à 120 Strips

Artikel-Nr.: 712155 | Biozym

89,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Farbe: frosted

Kategorie: 8er-Streifen

Produktlinie: Premium Line

Profil: High Profile

Verpackung: Beutel à 120 Strips

Volumen: 0.2 ml